

**Center for Veterinary Biologics
and
National Veterinary Services Laboratories
Testing Protocol**

**Supplemental Assay Method for Isolation of
Parainfluenza-3 Virus From Nasal Excretions**

Date: December 15, 2000
Supersedes: SAM 104 dated July 29, 1983
Number: MVSAM0104.01
Standard Requirement: 9 CFR, Part 113.309
Contact Person: Peg Patterson (515) 663-7334
Victor Becerra (515) 663-7468

Approvals:

Linn A. Wilbur, Head/Team Leader
Mammalian Virology Section
Date: _____

Ann L. Wieggers, Quality Assurance Manager
Date: _____

Randall L. Levings, Director
Center for Veterinary Biologics-Laboratory
Date: _____

United States Department of Agriculture
Animal and Plant Health Inspection Service
P. O. Box 844
Ames, IA 50010

Mention of trademark or proprietary product does not constitute a guarantee or warranty of the product by USDA and does not imply its approval to the exclusion of other products that may be suitable.

Supplemental Assay Method for Isolation of Parainfluenza-3 virus From Nasal
Excretions

Table of Contents

1. Introduction
 - 1.1 Background
 - 1.2 Keywords
2. Materials
 - 2.1 Equipment/instrumentation
 - 2.2 Reagents/supplies
3. Preparation for the test
 - 3.1 Personnel qualifications/training
 - 3.2 Preparation of equipment/instrumentation
 - 3.3 Preparation of reagents/control procedures
 - 3.4 Preparation of the sample
4. Performance of the test
5. Interpretation of the test results
6. Report of test results
7. References
8. Summary of revisions

Supplemental Assay Method for Isolation of Parainfluenza-3 virus From Nasal Excretions

1. Introduction

1.1 Background

This Supplemental Assay Method (SAM) describes an *in vitro* assay method using Madin Darby bovine kidney (MDBK) cells for the isolation of parainfluenza 3 virus (PI3V) from nasal excretions of calves. PI3V isolation is performed for immunogenicity studies on PI3 live virus vaccines in accordance with the Code of Federal Regulations, Title 9 (9 CFR).

1.2 Keywords

Parainfluenza 3 virus; PI3V; isolation, *in vitro*, CPE, HAd

2. Materials

2.1 Equipment/instrumentation

2.1.1 Incubator,¹ 36° ± 2°C, 5% ± 1% CO₂, high humidity

2.1.2 Water bath,² 36° ± 2°C

2.1.3 Pipettor,³ 200 µl and 1000 µl, and tips⁴

2.1.4 Vortex mixer⁵

2.1.5 Pipettor,⁶ 50-300 µl x 12 channel

2.1.6 Microscope,⁷ inverted light

¹ Model 3158, Forma Scientific, Inc., Box 649, Marietta, OH 45750-0649 or equivalent

² Model MW-1120A, Blue M Electric Co., 138th and Chatham St., Blue Island, IL or equivalent

³ Pipetman®, Rainin Instrument Co., Mack Rd., Box 4026, Woburn, MA 01888 or equivalent

⁴ Cat. No. YE-3R, Analytic Lab. Accessories, P.O. Box 345, Rockville Center, NY 11571 or equivalent

⁵ Vortex-3 Genie, Model G-560, Scientific Industries, Inc., Bohemia, NY 11716 or equivalent

⁶ Finnipette®, Labsystems OY, Pulttitie 9, 00810 Helsinki 81, Finland or equivalent

⁷ Model CK, Olympus America, Inc., 2 Corporate Center Dr., Melville, NY 11747-3157 or equivalent

Supplemental Assay Method for Isolation of Parainfluenza-3 virus From Nasal Excretions

2.1.7 Centrifuge⁸ and rotor⁹

2.1.8 Pipette-aid¹⁰

2.2 Reagent/supplies

2.2.1 PI3V Reference¹¹

2.2.2 MDBK cell line¹¹ found to be free of extraneous agents as tested by the 9 CFR.

2.2.3 Diluent Medium

2.2.3.1 9.61 g minimum essential medium (MEM) with Earle's salts without bicarbonate¹²

2.2.3.2 2.2 g sodium bicarbonate (NaHCO_3)¹³

2.2.3.3 Dissolve with 900 ml deionized water (DW).

2.2.3.4 Add 5.0 g lactalbumin hydrolysate or edamin¹⁴ to 10 ml DW. Heat to $60^\circ \pm 2^\circ\text{C}$ until dissolved and add to **Section 2.2.3.3** with constant mixing.

2.2.3.5 Q.S. to 1000 ml with DW; adjust pH to 6.8-6.9 with 2N hydrochloric acid (HCl).¹⁵

2.2.3.6 Sterilize through a 0.22- μm filter.¹⁶

⁸ Model J6-B, Beckman Coulter, P.O. Box 3100, Fullerton, CA 92834-3100 or equivalent

⁹ Type JS-4.0, Beckman Coulter or equivalent

¹⁰ Cat. No. 183, Drummond Scientific Co., 500 Pkwy., Broomall, PA 19008 or equivalent

¹¹ Reference quantities available upon request from the Center for Veterinary Biologics-Laboratory (CVB-L), P.O. Box 844, Ames, IA 50010 or equivalent

¹² Cat. No. 410-1500EF, Life Technologies, Inc., 8400 Helgeman Ct., Gaithersburg, MD 20884 or equivalent

¹³ Cat. No. S 5761, Sigma Chemical Co., P.O. Box 14508, St. Louis, MO 63178 or equivalent

¹⁴ Cat. No. 59102, Sheffield Products, P.O. Box 630, Norwick, NY 13815 or equivalent

¹⁵ Cat. No. 9535-01, J.T. Baker, Inc., 222 Red School Ln., Phillipsburg, NJ 08865 or equivalent

¹⁶ Cat. No. 12122, Gelman Sciences, 600 S. Wagner Rd., Ann Arbor, MI 48106 or equivalent

Supplemental Assay Method for Isolation of Parainfluenza-3 virus From Nasal Excretions

2.2.3.7 Aseptically add:

1. 10 ml L-glutamine¹⁷
2. 100 units/ml penicillin¹⁸
3. 50 µg/ml gentamicin sulfate¹⁹
4. 100 µg/ml streptomycin²⁰
5. 2.5 µg/ml amphotericin B²¹

2.2.3.8 Store at 4° ± 2°C.

2.2.4 Sample Transport Medium

2.2.4.1 200 ml Diluent Medium

2.2.4.2 60 units/ml penicillin

2.2.4.3 30 µg/ml gentamicin sulfate

2.2.4.4 60 µg/ml streptomycin

2.2.4.5 1.5 µg/ml amphotericin B

2.2.4.6 Store at 4° ± 2°C.

2.2.5 Growth Medium

2.2.5.1 900 ml of Diluent Medium

2.2.5.2 Aseptically add 100 ml gamma-irradiated fetal bovine serum (FBS).

2.2.5.3 Store at 4° ± 2°C.

2.2.6 Maintenance Medium

2.2.6.1 980 ml Diluent Medium

¹⁷ Cat. No. 320-503PE, Life Technologies, Inc., or equivalent

¹⁸ Cat. No. 0049-0530-28, Schering Laboratories, 2000-T Galloping Hill Rd., Kenilworth, NJ 07033 or equivalent

¹⁹ Cat. No. 0061-0464-04, Schering Laboratories or equivalent

²⁰ Cat. No. S-9137, Sigma Chemical Co. or equivalent

²¹ E.R. Squibb & Sons, Inc., 1 Squibb Dr., Cranberry, NJ 08512 or equivalent

Supplemental Assay Method for Isolation of Parainfluenza-3 virus From Nasal Excretions

2.2.6.2 Aseptically add 20 ml gamma-irradiated FBS.

2.2.6.3 Store at $4^{\circ} \pm 2^{\circ}\text{C}$.

2.2.7 Alsevers Solution

2.2.7.1 20.5 g dextrose ($\text{C}_6\text{H}_{12}\text{O}_6$)²²

2.2.7.2 8.0 g sodium citrate ($\text{Na}_3\text{C}_6\text{H}_5\text{O}_7 \cdot 2\text{H}_2\text{O}$)²³

2.2.7.3 4.2 g sodium chloride (NaCl)²⁴

2.2.7.4 0.55 g citric acid ($\text{C}_6\text{H}_8\text{O}_7 \cdot \text{H}_2\text{O}$)²⁵

2.2.7.5 Dissolve with 100 ml DW.

2.2.7.6 Filter through a 0.22- μm filter.

2.2.7.7 Store at $4^{\circ} \pm 2^{\circ}\text{C}$.

2.2.8 10X Phosphate buffered saline (10X PBS)

2.2.8.1 80.0 g sodium chloride

2.2.8.2 2.0 g potassium chloride (KCl)²⁶

2.2.8.3 2.0 g potassium phosphate, monobasic, anhydrous (KH_2PO_4)²⁷

2.2.8.4 Dissolve with 900 ml DW.

2.2.8.5 Add 11.5 g sodium phosphate dibasic, anhydrous (Na_2HPO_4)²⁸ to 50 ml DW, heat to

²²Cat. No. D12-500, Fisher Scientific Corp., 7711 Forbes Ave., Pittsburgh, PA 15219-4785 or equivalent

²³Cat. No. S279-500, Fisher Scientific Corp. or equivalent

²⁴Cat. No. 3624-01, J.T. Baker, Inc. or equivalent

²⁵Cat. No. A104-500, Fisher Scientific Corp. or equivalent

²⁶Cat. No. P217-500, Fisher Scientific Corp. or equivalent

²⁷Cat. No. 3246-01, J.T. Baker, Inc. or equivalent

²⁸Cat. No. 3828-01, J.T. Baker, Inc. or equivalent

Supplemental Assay Method for Isolation of Parainfluenza-3 virus From Nasal Excretions

60° ± 2°C until dissolved. Add to **Section 2.2.8.4** with constant mixing.

2.2.8.6 Q.S. to 1000 ml with DW.

2.2.8.7 Autoclave at 15 psi, 121° ± 2°C for 35 ± 5 min.

2.2.8.8 Store at 4° ± 2°C.

2.2.9 1X PBS (PBS)

2.2.9.1 100 ml 10X PBS

2.2.9.2 900 ml DW

2.2.9.3 Adjust the pH to 7.0-7.3 with 5N sodium hydroxide (NaOH).²⁹

2.2.9.4 Store at 4° ± 2°C.

2.2.10 Guinea pig red blood cells (RBCs) in an equal volume of Alsevers Solution

2.2.11 7.5 % Sodium Bicarbonate

2.2.11.1 7.5 g sodium bicarbonate

2.2.11.2 Q.S. to 100 ml DW.

2.2.11.3 Autoclave at 15 psi, 121° ± 2°C for 30 ± 5 min.

2.2.11.4 Store at 4° ± 2°C.

2.2.12 Trypsin Versene (TV)

2.2.12.1 8.0 g sodium chloride

²⁹ Cat. No. SS256-500, Fisher Scientific Co. or equivalent

Supplemental Assay Method for Isolation of Parainfluenza-3 virus From Nasal Excretions

- 2.2.12.2 0.40 g potassium chloride
- 2.2.12.3 0.58 g sodium bicarbonate
- 2.2.12.4 0.50 g irradiated trypsin³⁰
- 2.2.12.5 0.20 g versene or disodium salt ethylenediaminetetraacetic acid (EDTA)³¹
- 2.2.12.6 1.0 g dextrose
- 2.2.12.7 0.4 ml 0.5% phenol red³²
- 2.2.12.8 Dissolve with 1000 ml DW.
- 2.2.12.9 Adjust pH to 7.3 with 7.5% Sodium Bicarbonate.
- 2.2.12.10 Filter through a 0.22- μ m filter.
- 2.2.12.11 Aseptically dispense into 100 ml amounts and store at $-20^{\circ} \pm 4^{\circ}\text{C}$.
- 2.2.13 Tissue culture plates,³³ 96-well
- 2.2.14 Polystyrene tubes, 12 x 75 mm³⁴ and 17 x 100³⁵
- 2.2.15 Conical tube,³⁶ 50 ml
- 2.2.16 Serological pipettes,³⁷ 10 ml
- 2.2.17 Cotton-tipped applicators,³⁸ sterile

³⁰ Cat. No. 0152-15-9, DIFCO Laboratories, P.O. Box 331058, Detroit, MI 48232-0758

³¹ Cat. No. E 5134, Sigma Chemical Co. or equivalent

³² Cat. No. P0290, Sigma Chemical Co. or equivalent

³³ Cat. No. 3596, Costar Corning Inc., 1 Alewife Center, Cambridge, MA 02140 or equivalent

³⁴ Falcon® 2058, Becton Dickinson Labware or equivalent

³⁵ Falcon® 2057, Becton Dickinson Labware or equivalent

³⁶ Cat. No. 62.547, Sarstedt, Inc. P.O. Box 468, Newton, NC 28658-0468 or equivalent

³⁷ Falcon® 7530, Becton Dickinson Labware or equivalent

³⁸ Cat. No. 22-9588, The CITMED Corp., 18601 S. Main St., Citronelle, AL 36522 or equivalent

Supplemental Assay Method for Isolation of Parainfluenza-3 virus From Nasal Excretions

2.2.18 Reagent reservoir³⁹

3. Preparation for the test

3.1 Personnel qualifications/training

Personnel must have training in cell culture techniques, virus isolation, the principles of aseptic technique, and animal care and handling techniques.

3.2 Preparation of equipment/instrumentation

On the day of test initiation, set a water bath at $36^{\circ} \pm 2^{\circ}\text{C}$.

3.3 Preparation of reagents/control procedures

3.3.1 Sample Transport Tube. Prior to Nasal Excretion Sample collection, fill an appropriate number of 12 x 75-mm polystyrene tubes with 3 ml of sample transport medium for each calf on test. A minimum of 1 tube per calf per day of anticipated collection is prepared.

3.3.2 Nasal Excretion Sample collections

Take Nasal Excretion Samples from each calf by inserting a sterile cotton-tipped applicator several inches into each nasal passage. Immediately immerse the two applicators laden with nasal excretion material into a Sample Transport Tube; labeled appropriately with calf identifier and date of collection. Freeze and store the Sample Transport Tube at $-70^{\circ} \pm 5^{\circ}\text{C}$ until virus isolation test is conducted.

3.3.3 Preparation of MDBK cell culture plate (Test Plate)

³⁹Cat. No. 4870, Costar Corning Inc., 1 Riverfront Plaza, Corning, NY 14831 or equivalent

Supplemental Assay Method for Isolation of Parainfluenza-3 virus From Nasal Excretions

Cells are prepared from healthy, confluent MDBK cells, that are maintained by passing weekly. Two days prior to test initiation and 2 days prior to the 2nd passage, seed 200 µl/well of $10^{5.4}$ to $10^{5.6}$ cells/ml, suspended in Growth Medium into all wells of a 96-well cell culture plate with a multichannel pipettor. At least 2 Nasal Excretion Samples are tested on 1 plate. Cells are removed from the growth containers by using TV solution. The Test Plate is incubated at $36^{\circ} \pm 2^{\circ}\text{C}$ in a CO_2 incubator for 36 ± 12 hours. Cells should be 80% confluent prior to use.

3.3.4 Preparation of PI3V Reference Control

3.3.4.1 On the day of test initiation, rapidly thaw a vial of PI3V Reference in a $36^{\circ} \pm 2^{\circ}\text{C}$ water bath and make tenfold dilutions.

1. Using a 10-ml serological pipette, dispense 4.5 ml of Dilution Medium in an appropriate number of 17 x 100-mm polystyrene tubes to bracket the expected titer endpoint according to the CVB-L Reference and Reagent sheet. Label tubes (for example: 7 tubes, labeled 10^{-1} to 10^{-7} respectively).
2. Using a 500 µl pipettor, transfer 500 µl of PI3V Reference to the first tube labeled 10^{-1} ; mix by vortexing.
3. Using a new pipette tip, transfer 500 µl from the 10^{-1} labeled tube (**Section 3.3.4.1.2**) to the 10^{-2} tube; mix by vortexing.
4. Repeat **Section 3.3.4.1.3** for each of the subsequent dilutions, transferring 500 µl of the previous dilution to the next dilution tube, until the tenfold dilution series is completed.

Supplemental Assay Method for Isolation of Parainfluenza-3 virus From Nasal Excretions

3.3.5 Preparation of 0.5% RBCs Suspension for the hemadsorption (HAd) test

3.3.5.1. Transfer RBCs to a 50-ml conical tube.

3.3.5.2. Q.S. to 50 ml with Alsevers Solution; mix by inverting several times.

3.3.5.3. Centrifuge at $400 \times g$ for 15 ± 5 min (1500 rpm in a J6B centrifuge with a JS-4.0 rotor).

3.3.5.4. Remove supernatant and buffy coat by aspirating with a 10-ml serological pipette.

3.3.5.5. Repeat **Section 3.3.5.2 through Section 3.3.5.3** for a total of 3 washes, removing the supernatant each time.

3.3.5.6. Make a 0.5% RBC suspension by pipetting 500 μ l of packed RBCs into 100 ml of PBS; mix by inverting.

3.3.5.7. Store at $4^\circ \pm 2^\circ\text{C}$; use within a wk of collection of the RBCs.

3.4 Preparation of the sample

3.4.1 On the day of the test initiation, thaw the Sample Transport Tubes containing the nasal swabs, express the fluid from the cotton by rotating and pressing the applicators against the tube wall, and vortex the media in the tubes to release virus from the swab. Discard the used cotton applicators, after expressing any fluid.

3.4.2 Centrifuge the sample at $500 \times g$ (2000 rpm in a J6B centrifuge with a JS-4.0 rotor) for 20 ± 5 min. Transfer 2 ml of the supernatant from each sample into a new, labeled 12 x 75-mm polystyrene tube.

Supplemental Assay Method for Isolation of Parainfluenza-3 virus From Nasal Excretions

3.4.3 Keep the supernatant samples on ice until time of inoculation.

3.4.4 Store the remaining supernatant sample at $-70^{\circ} \pm 5^{\circ}\text{C}$.

4. Performance of the test

4.1 Label the Test Plate and aseptically decant the Growth Medium into a suitable container.

4.2 Inoculate 5 wells with 25 μl /well of each supernatant sample. Change tips between each sample.

4.3 Inoculate 5 wells/dilution with 25 μl /well of the PI3V Reference Control (with dilutions 10^{-7} through 10^{-4} from the example in **Section 3.3.4.1**). Tip changes are not necessary if pipetting from the most dilute (10^{-7}) to the most concentrated (10^{-4}).

4.4 Maintain 5 wells as uninoculated cell culture controls.

4.5 Incubate the Test Plate at $36^{\circ} \pm 2^{\circ}\text{C}$ in a CO_2 incubator for 60 ± 10 min for virus adsorption.

4.6 Add 200 μl /well of Maintenance Medium to the Test Plate with a multichannel pipettor.

4.7 Incubate the Test Plate at $36^{\circ} \pm 2^{\circ}\text{C}$ in a CO_2 incubator for 96 ± 24 hr. Observe daily for PI3V CPE which is characterized by cell fusion.

4.8 Harvest samples from each well in which typical PI3V infection is observed into labeled 12 x 75-mm polystyrene tubes. Samples may be harvested at time of observation. Freeze at $-70^{\circ} \pm 5^{\circ}\text{C}$ until next passage.

4.9 At the end of incubation, fluids from wells showing no CPE are passaged onto fresh Test Plates prepared according to **Section 3.3.3**, by repeating **Sections 4.2 through 4.7**. In

Supplemental Assay Method for Isolation of Parainfluenza-3 virus From Nasal Excretions

addition, the harvested samples that displayed PI3V infection are thawed and passed. A PI3V Reference Control prepared as described in **Sections 3.3.4 and 4.3**, respectively, is included. At the end of the 96 ± 24 hr incubation, an HAd test is performed on all second passage wells to confirm PI3V infection.

4.10 HAd Test

4.10.1 Decant Maintenance Medium from the inoculated Test Plates (**Section 4.9**) in a suitable container.

4.10.2 Rinse all cells of the Test Plates by filling with PBS; decant.

4.10.3 Dispense 200 μ l/well of a 0.5% of the RBC Suspension into all test wells.

4.10.4 Incubate the Test Plates 15 ± 5 min at room temperature ($23^{\circ} \pm 2^{\circ}\text{C}$).

4.10.5 Decant the RBC Suspension from the Test Plates and repeat **Section 4.10.2** for a total of 2 washes.

4.10.6 Examine the Test Plate monolayers on an inverted light microscope at 100X magnification, and record the results. Wells containing one or more RBC clusters adhering to the cell monolayer are considered to be positive for PI3V.

4.11 Calculate the PI3V endpoint of the PI3V Reference Control using the method of Spearman-Kärber as commonly modified. The titers are expressed as \log_{10} 50% tissue culture infective dose (TCID_{50}) of the test wells.

Example:

10^{-4} dilution of PI3V Reference = 5/5 wells CPE/HAd +
 10^{-5} dilution of PI3V Reference = 5/5 wells CPE/HAd +
 10^{-6} dilution of PI3V Reference = 2/5 wells CPE/HAd +
 10^{-7} dilution of PI3V Reference = 0/5 wells CPE/HAd +

Supplemental Assay Method for Isolation of Parainfluenza-3 virus From Nasal Excretions

Titer = $(X - d/2 + [d * S])$ where:
X = \log_{10} of lowest dilution (4)
d = \log_{10} of dilution factor (1)
S = sum of proportion of CPE/HAd positive

$$\frac{(5 + 5 + 2)}{5} = \frac{12}{5} = 2.4$$

PI3V Reference Control titer = $(4 - 1/2 + (1 * 2.4)) = 5.9$

Titer of the PI3V Reference Control is $10^{5.9}$ TCID₅₀.

5. Interpretation of the test results

5.1 The uninoculated cell controls should not exhibit degradation or cloudy media that would indicate contamination.

5.2 The calculated HAd₅₀ titer of the PI3V Reference Control must fall within plus or minus 2 standard deviations (± 2 SD) of its mean titer, as established from a minimum of 10 previously determined titers.

5.3 Wells of Nasal Excretion Samples that tested HAd positive (regardless if CPE developed) are considered positive for PI3V, whereas those wells of samples that did not develop CPE and tested HAd negative after the second serial passage in MDBK cells are considered to be negative for PI3V.

5.4 Any Nasal Excretion Sample with 1 or more positive wells is considered to be positive for PI3V. A Nasal Excretion Sample is only considered negative if all wells inoculated from the Nasal Excretion Sample are negative for PI3V.

Supplemental Assay Method for Isolation of Parainfluenza-3 virus From Nasal Excretions

6. Report of test results

6.1 Results for each Nasal Excretion Sample are reported as positive or negative PI3V isolation.

7. References

7.1 Code of Federal Regulations, Title 9, Part 113.309, U.S. Government Printing Office, Washington, D.C., 2000.

7.2 Cottral, G.E., (Ed.). *Manual of standardized methods for veterinary microbiology*. Comstock Publishing Associates, Ithaca and London, 1978, pg.731.

7.3 Finney, D.J. *Statistical method in biological assay*. Griffin, London. 3rd ed., 1978, pg. 508.

8. Summary of revisions

This document was rewritten to meet the current NVSL/CVB QA requirements, to clarify practices currently in use in the CVB-L, and to provide additional detail. Changes made from the previous protocol are as follows:

- 1) use of Madin-Darby bovine kidney-A cells instead of primary embryonic kidney cell cultures,
- 2) replacement of the 24-well cell culture plates by 96-well plates,
- 3) use of minimum essential medium from brain heart infusion broth for transport medium,
- 4) a change from 5 days to 4 ± 1 day incubation.