

**Center for Veterinary Biologics  
and  
National Veterinary Services Laboratories  
Testing Protocol**

**Supplemental Assay Method for Titration of  
Pseudorabies Virus in Vaccines**

Date: October 29, 1998

Supersedes: SAM 118, Dated: March 1, 1989

Number: MVSAM0118.01

Standard Requirement: 9 CFR 113.318

Contact Person: Peg A. Patterson, (515) 239-8334  
Larry R. Ludemann, (515) 239-8264

Approvals:

\_\_\_\_\_  
Linn A. Wilbur, Head/Team Leader  
Mammalian Virology Section

Date: \_\_\_\_\_

\_\_\_\_\_  
Richard E. Hill, Quality Assurance Manager

Date: \_\_\_\_\_

\_\_\_\_\_  
Randall L. Levings, Director  
Center for Veterinary Biologics-Laboratory

Date: \_\_\_\_\_

United States Department of Agriculture  
Animal and Plant Health Inspection Service  
P. O. Box 844  
Ames, IA 50010

Mention of trademark or proprietary product does not constitute a guarantee or warranty of the product by USDA and does not imply its approval to the exclusion of other products that may be suitable.

Supplemental Assay Method for Titration of Pseudorabies Virus in Vaccines

Table of Contents

1. Introduction
2. Materials
  - 2.1 Equipment/instrumentation
  - 2.2 Reagents/supplies
3. Preparation for the test
  - 3.1 Personnel qualifications/training
  - 3.2 Preparation of equipment/instrumentation
  - 3.3 Preparation of reagents/control procedures
  - 3.4 Preparation of the sample
4. Performance of the test
5. Interpretation of the test results
6. Report of test results
7. References
8. Changes

Supplemental Assay Method for Titration of Pseudorabies Virus in Vaccines

## 1. Introduction

This is an *in vitro* test method which utilizes viral plaque forming units (PFU) in a cell culture system to titer pseudorabies virus (PRV) in modified-live vaccines.

## 2. Materials

### 2.1 Equipment/instrumentation

2.1.1 36° ± 2°C, 5% ± 1% CO<sub>2</sub>, high humidity incubator<sup>1</sup> meeting the requirements of the current version of GDOCSOP004

2.1.2 Vortex mixer<sup>2</sup>

2.1.3 Blender<sup>3</sup>

2.1.4 Micropipettors: 200 µl and 1000 µl single channel<sup>4</sup>

2.1.6 Water bath<sup>5</sup>

2.1.7 Self-refilling, 2-ml repetitive syringe<sup>6</sup>

### 2.2 Reagents/supplies

2.2.1 PRV Reference Virus,<sup>7</sup> Shope strain

2.2.2 Madin-Darby bovine kidney<sup>8</sup> (MDBK) cells

---

<sup>1</sup> Model 3158, Forma Scientific Inc., Box 649, Marietta, OH 45750-0649 or equivalent

<sup>2</sup> Vortex-2 Genie, Model G-560, Scientific Industries Inc., 70 Orville Dr., Bohemia, NY 11716 or equivalent

<sup>3</sup> Waring blender, Cat. No. 14-509-35, Fisher Scientific Inc., 319 W. Ontario, Chicago, IL 60610 or equivalent

<sup>4</sup> Pipetman, Rainin Instrument Co., Mack Rd., Box 4026, Woburn, MA 01888 or equivalent

<sup>5</sup> Cat. No. 15-461-10, Fisher Scientific Inc. or equivalent

<sup>6</sup> Wheaton, Cat. No. 13-689-50C, Fisher Scientific Inc. or equivalent

<sup>7</sup> Reference quantities are available upon request from the Center for Veterinary Biologics-Laboratory (CVB-L), P.O. Box 844, Ames, IA 50010 or equivalent

<sup>8</sup> Cat. No. ATCC CCL-22, American Type Culture Collection, 12301 Parklawn Dr., Rockville, MD 20852-1776

Supplemental Assay Method for Titration of Pseudorabies Virus in Vaccines

2.2.3 Minimum Essential Medium (MEM)

2.2.3.1 9.61 g MEM<sup>9</sup>

2.2.3.2 2.2 g sodium bicarbonate<sup>10</sup> (NaHCO<sub>3</sub>)

2.2.3.3 Q.S. to 1000 ml with deionized water (DW), adjust pH to 6.8-6.9 with 2N hydrochloric acid (HCl).<sup>11</sup>

2.2.3.4 Sterilize through 0.22-µm filter.<sup>12</sup>

2.2.3.5 Aseptically add:

1. 10 ml L-glutamine<sup>13</sup>
2. 5 ml lactalbumin hydrolysate or edamin<sup>14</sup>
3. 100 units/ml penicillin<sup>15</sup>
4. 50 µg/ml gentamicin sulfate<sup>16</sup>
5. 100 µg/ml streptomycin<sup>17</sup>

2.2.3.6 Store at 4° ± 2°C.

2.2.4 Growth Medium

2.2.4.1 900 ml of MEM

---

<sup>9</sup> MEM with Earle's salts without sodium bicarbonate, Cat. No. 410-1500EF, Life Technologies, Inc., 8400 Helgerman Ct., Gaithersburg, MD 20884 or equivalent

<sup>10</sup> Cat. No. S 5761, Sigma Chemical, Inc., P.O. Box 14508, St. Louis, MO 63178 or equivalent

<sup>11</sup> Cat. No. 9535-01, J.T. Baker, Inc., 222 Red School Ln., Phillipsburg, NJ 08865 or equivalent

<sup>12</sup> Cat. No. 12122, Gelman Sciences, 600 S. Wagner Rd., Ann Arbor, MI 48106 or equivalent

<sup>13</sup> L-glutamine-200 mM (100X), liquid, Cat. No. 320-503PE, Life Technologies, Inc. or equivalent

<sup>14</sup> Edamin S, Cat. No. 59102, Sheffield Products, P.O. Box 630, Norwick, NY 13815 or equivalent

<sup>15</sup> Cat. No. 0049-0530-28, Schering Laboratories, 2000-T Galloping Hill Rd., Kenilworth, NJ 07033 or equivalent

<sup>16</sup> Cat. No. 0061-0464-04, Schering Laboratories or equivalent

<sup>17</sup> Cat. No. S-9137, Sigma Chemical, Inc. or equivalent

Supplemental Assay Method for Titration of Pseudorabies Virus in Vaccines

2.2.4.2 Aseptically add 100 ml fetal bovine serum (FBS), heat-inactivated at  $56^{\circ} \pm 2^{\circ}\text{C}$  for  $30 \pm 5$  min.

2.2.4.3 Store at  $4^{\circ} \pm 2^{\circ}\text{C}$ .

2.2.5 2X Medium

2.2.5.1 100 ml 10X MEM<sup>18</sup>

2.2.5.2 2.2 g sodium bicarbonate

2.2.5.3 340 ml DW

2.2.5.4 Sterilize through 0.22- $\mu\text{m}$  filter.

2.2.5.5 Aseptically add:

1. 100 units/ml penicillin
2. 50  $\mu\text{g}/\text{ml}$  gentamicin sulfate
3. 100  $\mu\text{g}/\text{ml}$  streptomycin
4. 50 ml FBS

2.2.5.6 Store at  $4^{\circ} \pm 2^{\circ}\text{C}$ .

2.2.6 2% Tragacanth Gum (Trag)

2.2.6.1 20 g Trag<sup>19</sup>

2.2.6.2 1000 ml DW

2.2.6.3 Mix vigorously small amounts at a time with a blender set on high.

2.2.6.4 Pour 500 ml each into 1000-ml media bottles.

2.2.6.5 Sterilize by autoclaving at 15 psi for 30 min.

---

<sup>18</sup> MEM with Earle's salts 10X liquid, Cat. No. 410-11430, Life Technologies, Inc. or equivalent

<sup>19</sup> Acros AC42138-5000, Fisher Scientific, Inc. or equivalent

Supplemental Assay Method for Titration of Pseudorabies Virus in Vaccines

2.2.6.6 Store at  $4^{\circ} \pm 2^{\circ}\text{C}$ .

2.2.7 Overlay Medium

2.2.7.1 Mix equal volumes of 2X Medium and 2% Trag.

2.2.7.2 Store at  $4^{\circ} \pm 2^{\circ}\text{C}$ .

2.2.8 70% Ethyl Alcohol

2.2.8.1 74 ml ethyl alcohol<sup>20</sup>

2.2.8.2 26 ml DW

2.2.8.3 Store at room temperature (RT) ( $23^{\circ} \pm 2^{\circ}\text{C}$ ).

2.2.9 Crystal Violet Stain

2.2.9.1 7.5 g crystal violet<sup>21</sup>

2.2.9.2 50 ml 70% ethyl alcohol

2.2.9.3 Dissolve crystal violet in alcohol.

2.2.9.4 250 ml formaldehyde<sup>22</sup>

2.2.9.5 Q.S. to 1000 ml with DW, filter through filter paper.<sup>23</sup>

2.2.9.6 Store at RT.

2.2.10 Tissue culture plates, 6 well<sup>24</sup>

2.2.11 12x75-mm polystyrene tubes<sup>25</sup>

---

<sup>20</sup> Denatured, 190 proof, Cat. No. 7018, J.T. Baker, Inc. or equivalent

<sup>21</sup> Cat. No. C0775, Sigma Chemical Co. or equivalent

<sup>22</sup> Cat. No. F79, Fisher Scientific, Inc. or equivalent

<sup>23</sup> Whatman #1, Cat. No. 1001, Fisher Scientific, Inc. or equivalent

<sup>24</sup> Falcon 3046, Becton Dickinson Labware, Becton Dickinson & Co., 2 Oak Park, Bedford, MA 01730 or equivalent

<sup>25</sup> Falcon 2058, Becton Dickinson Labware or equivalent

**Supplemental Assay Method for Titration of Pseudorabies Virus in Vaccines**

**2.2.12 10-ml syringes<sup>26</sup> and needles<sup>27</sup>**

**3. Preparation for the test**

**3.1 Personnel qualifications/training**

Personnel must have training in the immunologic basis of virus titration assays, cell culture maintenance, and in the principles of aseptic techniques.

**3.2 Preparation of equipment/instrumentation**

Set the water bath at  $36^{\circ} \pm 2^{\circ}\text{C}$ .

**3.3 Preparation of reagents/controls**

**3.3.1. Two days prior to test performance**

Seed 6-well tissue culture plates with MDBK cells, in Growth Medium, at a cell count that will produce a monolayer after  $48 \pm 6$  hr of incubation at  $36^{\circ} \pm 2^{\circ}\text{C}$ . These become the MDBK plates. Growth Medium is changed if excess acidity of the medium is observed as indicated by a change from red to yellow of Growth Medium.

**3.3.2 On day of test performance**

Rapidly thaw a vial of PRV Reference Virus in a water bath at  $36^{\circ} \pm 2^{\circ}\text{C}$ . Dilute the virus in 4.5 ml MEM to contain 15-40 PFU/100  $\mu\text{l}$ .

**3.4 Preparation of the sample (on day of test performance)**

**3.4.1** Rehydrate a vial of the Test Vaccine according to the manufacturer's instructions with a 10-ml syringe and 20-gauge needle. Allow to sit for  $15 \pm 5$  min at RT.

---

<sup>26</sup> Luer-Lok®, Cat. No. 309604, Becton Dickinson Labware or equivalent

<sup>27</sup> 20 gauge, 1.5 in, Cat. No. 250107, Becton Dickinson Labware or equivalent

**Supplemental Assay Method for Titration of Pseudorabies Virus in Vaccines**

**3.4.2** Prepare serial 10-fold dilutions of test vaccine. Serial 10-fold dilutions may be made as follows:

**3.4.2.1** Place 4.5 ml of MEM with the repetitive syringe into labeled 12x75-mm polystyrene tubes.

**3.4.2.2** Pipet 500 µl of test vaccine to the 10<sup>-1</sup> tube, mix by vortexing. Discard pipet tip.

**3.4.2.3** Repeat step **3.4.2.2** to the remaining tubes. Continue as needed (10<sup>-2</sup> to 10<sup>-5</sup>), transferring 500 µl from previous tube to the next dilution.



Supplemental Assay Method for Titration of Pseudorabies Virus in Vaccines

**4. Performance of the test**

- 4.1 Decant the Growth Media from MDBK plates.
- 4.2 Pipet 100  $\mu$ l/well from each dilution of test vaccine to 2 wells of a MDBK plate. Mix by gentle swirling.
- 4.3 Pipet 100  $\mu$ l/well of the diluted Reference Virus Control to 2 wells of 1 MDBK plate. Mix by gentle swirling.
- 4.4 Maintain 2 or more wells as uninoculated cell culture controls.
- 4.5 Incubate inoculated plates at  $36^{\circ} \pm 2^{\circ}\text{C}$  in a  $\text{CO}_2$  atmosphere for  $60 \pm 10$  min for virus adsorption.
- 4.6 Add 3.0 ml/well of Overlay Medium (see section 2.2.8) to the plates. Discard any unused, warmed Overlay Medium.
- 4.7 Incubate the MDBK plates undisturbed at  $36^{\circ} \pm 2^{\circ}\text{C}$  in a  $\text{CO}_2$  atmosphere for  $96 \pm 12$  hr.
- 4.8 At the end of incubation, decant Overlay Medium. Pipet 2 ml of the Crystal Violet Stain (see section 2.2.10) into each well of the plates using the repetitive syringe.
- 4.9 Allow plates to stand at RT for  $15 \pm 5$  min.
- 4.10 Discard the Crystal Violet Stain down a sink. Wash the cell monolayers by dipping each plate several times in a container of running water from the cold faucet. Allow to air dry.
- 4.11 PFU counting
  - 4.11.1 The PFU are visible as clear, circular areas in the cell monolayer where the cells have been destroyed by the virus.

**Supplemental Assay Method for Titration of Pseudorabies Virus in Vaccines**

**4.11.2** Count the number of PFU for each well.

**4.11.2.1** Average the number of PFU between the duplicate wells for each test vaccine dilution.

**4.11.2.2** Average the number of PFU between the 2 wells of the Reference Virus Control wells.

**4.11.3** Determine the virus titers and express as PFU/dose. Take the dilution of a test vaccine that contains an average of at least 30 PFU's. Example:

Log <sub>10</sub> of plaque count (30)	1.48
Log <sub>10</sub> of dilution counted (10 <sup>-3</sup> )	3.00
Log <sub>10</sub> of 2-ml dose factor (20)	1.30
<hr/>	
Virus titer/dose (total)	5.78

The test vaccine contains 10<sup>5.78</sup> PFU/2-ml dose.

**5. Interpretation of the test results**

**5.1** The test is invalid if visible contamination is observed in all dilutions of a Test Vaccine.

**5.2** For a valid assay, the Reference Virus Control must have an average PFU count between 15-40.

**5.3** Any test not meeting the criteria of **5.1** and **5.2** is considered a NO TEST and may be repeated without prejudice.

**5.4** One plaque represents a single infective unit whereas the 50% endpoint infective dose (one ID<sub>50</sub>) is statistically equivalent to a theoretical 0.69 of an infective unit. Fifty percent endpoints will be 1.44 times those expressed as PFU/unit of inoculation. Therefore, to express PFU titer as 50% tissue culture infective dose (TCID<sub>50</sub>) titer, multiply the PFU by 1.44 or add 0.16 to the log<sub>10</sub> value of the PFU titer. In the example above, the TCID<sub>50</sub> would be 5.78 + 0.16 = 5.94 or 10<sup>5.94</sup> TCID<sub>50</sub>/2-ml dose.

**5.5** If the validity requirements are not met, then the assay is considered a NO TEST and can be retested without prejudice.

**Supplemental Assay Method for Titration of Pseudorabies Virus in Vaccines**

**5.6** If the validity requirements are met and the titer of the vaccine is greater than or equal to the titer contained in the filed outline of production for the product under test, the product is considered satisfactory.

**5.7** If the validity requirements are met but the titer of the Test Vaccine is lower than the required minimum, it must be retested according to 9 CFR 113.8.

**6. Report of test results**

**6.1** Test results are reported as the virus titer in TCID<sub>50</sub>/dose.

**6.2** Record all test results on the test record.

**7. References**

**7.1** Conrath TB.: Handbook of Microtiter Procedures. 1972. Clinical and Research Applications Laboratory, Cooke Engineering Company, Alexandria, VA.

**8. Changes**

**8.1** This document was rewritten to meet the current NVSL/CVB QA requirements, to clarify practices currently in use in the CVB-L, and to provide additional detail. No significant changes were made from the previous protocol.