

**Center for Veterinary Biologics
and
National Veterinary Services Laboratories
Testing Protocol**

**Supplemental Assay Method for the Titration of Canine
Parvovirus in Cell Culture**

Date: December 22, 1999

Supersedes: August 1, 1987

Number: MVSAM0316.01

Standard Requirement: 9 CFR, Part 113.317

Contact Person: Larry R. Ludemann (515) 663-7264
Marsha J. Hegland (515) 663-7659

Approvals:

Linn A. Wilbur, Head/Team Leader
Mammalian Virology Section

Date: _____

Ann L. Wieggers, Quality Assurance Manager

Date: _____

/s/ Randall Levings _____ Date: 12/22/99
Randall L. Levings, Director
Center for Veterinary Biologics-Laboratory

United States Department of Agriculture
Animal and Plant Health Inspection Service
P. O. Box 844
Ames, IA 50010

Mention of trademark or proprietary product does not constitute a guarantee or warranty of the product by USDA and does not imply its approval to the exclusion of other products that may be suitable.

Supplemental Assay Method for the Titration of Canine Parvovirus
in Cell Culture

Table of Contents

1. Introduction
 - 1.1 Background
 - 1.2 Keywords
2. Materials
 - 2.1 Equipment/instrumentation
 - 2.2 Reagents/supplies
3. Preparation for the test
 - 3.1 Personnel qualifications/training
 - 3.2 Preparation of equipment/instrumentation
 - 3.3 Preparation of reagents/control procedures
 - 3.4 Preparation of the sample
4. Performance of the test
5. Interpretation of the test results
6. Report of test results
7. References
8. Summary of revisions

Supplemental Assay Method for the Titration of Canine Parvovirus
in Cell Culture

1. Introduction

1.1 Background

This Supplemental Assay Method (SAM) is an *in vitro* test method for assaying modified-live canine parvovirus (CPV) vaccines for viral content. The method uses the Crandall feline kidney (CRFK) cell line as the test system. Presence or absence of CPV is determined by staining inoculated cell cultures by an indirect fluorescent antibody (IFA) technique.

1.2 Keywords

Canine parvovirus; CPV; IFA; potency test; TCID₅₀; viral titration; *in vitro*

2. Materials

2.1 Equipment/instrumentation

2.1.1 Incubator,¹ 36° ± 2°C, high humidity, 5 ± 1% CO₂ (CO₂ incubator) meeting the requirements of the current version of GDOCSOP004

2.1.2 Incubator,² 36° ± 2°C, aerobic

2.1.3 Water bath,³ 36° ± 2°C

2.1.4 Microscope,⁴ ultraviolet (UV) light

2.1.5 Vortex mixer⁵

2.1.6 Micropipettor,⁶ 200 µl with tips⁷

¹ Model 3336, Forma Scientific, Inc., P.O. Box 649, Marietta, OH 45750 or equivalent

² Model 2, Precision Scientific, 3737 West Cortland St., Chicago, IL 60647 or equivalent

³ Cat. No. 66648, Precision Scientific or equivalent

⁴ Model BH2, Olympus America, Inc., 2 Corporate Ctr. Dr., Melville, NY 11747 or equivalent

⁵ Vortex-2 Genie, Model G-560, Scientific Industries, Inc., 70 Orville Dr., Bohemia, NY 11716 or equivalent

⁶ Cat. No. P-200, Rainin Instrument Co., P.O. Box 4026, Mack Rd., Woburn, MA 01801-4628 or equivalent

⁷ Cat. No. YE-3R, Analytic Lab Accessories, P.O. Box 345, Rockville Centre, NY 11571 or equivalent

Supplemental Assay Method for the Titration of Canine Parvovirus
in Cell Culture

2.1.7 Microscope slide glass staining dish with rack
(glass staining dish)⁸

2.2 Reagents/supplies

2.2.1 CPV Reference,⁹ KB5 strain

2.2.2 CRFK cell culture,¹⁰ free of extraneous agents
as tested by the Code of Federal Regulations,
Title 9 (9 CFR)

2.2.3 Feline panleukopenia virus antiserum¹¹ (FPV
Antiserum)

2.2.4 Minimum essential medium (MEM)

2.2.4.1 9.61 g MEM with Earle's salts without
bicarbonate¹²

2.2.4.2 2.2 g sodium bicarbonate (NaHCO₃)¹³

2.2.4.3 Dissolve with 900 ml deionized water
(DW).

2.2.4.4 Add 5.0 g lactalbumin hydrolysate or
edamine¹⁴ to 10 ml DW. Heat to 60° ± 2°C until
dissolved. Add to **Section 2.2.5.3** with constant
mixing.

2.2.4.5 Q.S. to 1000 ml with DW; adjust pH to
6.8-6.9 with 2N hydrochloric acid (HCl).¹⁵

2.2.4.6 Sterilize through a 0.22-µm filter.¹⁶

⁸ Cat. No. 121, Shandon Lipshaw, 171 Industry Dr., Pittsburgh, PA 15275 or equivalent

⁹ Seed quantities available upon request from the Center for Veterinary Biologics-Laboratory
(CVB-L), P.O. Box 844, Ames, IA 50010 or equivalent

¹⁰ CCL-94, American Type Culture Collection, 12301 Parklawn Dr., Rockville, MD 20852-1776

¹¹ Reference quantities available upon request from the CVB-L or equivalent

¹² Cat. No. 410-1500EF, Life Technologies, Inc., 8400 Helgeman Ct., Gaithersburg, MD 20884 or
equivalent

¹³ Cat. No. S-5761, Sigma Chemical Co., P.O. Box 14508, St. Louis, MO 63178 or equivalent

¹⁴ Edamine, Cat. No. 59102, Sheffield Products, P.O. Box 630, Norwick, NY 13815 or equivalent

¹⁵ Cat. No. 9535-01, J.T. Baker, Inc., 222 Red School Ln., Phillipsburg, NJ 08865 or equivalent

¹⁶ Cat. No. 12122, Gelman Sciences, 600 S. Wagner Rd., Ann Arbor, MI 48106 or equivalent

Supplemental Assay Method for the Titration of Canine Parvovirus
in Cell Culture

2.2.4.7 Aseptically add:

1. 25 units/ml penicillin¹⁷
2. 50 µg/ml gentamicin sulfate¹⁸
3. 100 µg/ml streptomycin¹⁹

2.2.4.8 Store at 4° ± 2°C.

2.2.5 Growth Medium

2.2.5.1 920 ml of MEM

2.2.5.2 Aseptically add:

1. 70 ml gamma irradiated fetal bovine serum (FBS)
2. 10 ml L-glutamine²⁰

2.2.6 Dulbecco's phosphate buffered saline (DPBS)

2.2.6.1 8.0 g sodium chloride (NaCl)²¹

2.2.6.2 0.2 g potassium chloride (KCl)²²

2.2.6.3 0.2 g potassium phosphate, monobasic, anhydrous (KH₂PO₄)²³

2.2.6.4 0.1 g magnesium chloride, hexahydrate (MgCl₂•6H₂O)²⁴

2.2.6.5 Dissolve reagents with 900 ml DW.

¹⁷Cat. No. 0049-0530-28, Schering Laboratories, 2000-T Galloping Hill Rd., Kenilworth, NJ 07033 or equivalent

¹⁸Cat. No. 0061-0464-04, Schering Laboratories or equivalent

¹⁹Cat. No. S-9137, Sigma Chemical Co. or equivalent

²⁰L-glutamine-200 mM (100X), liquid, Cat. No. 320-503PE, Life Technologies, Inc. or equivalent

²¹Cat. No. 3624-01, J.T. Baker, Inc. or equivalent

²²Cat. No. P217-500, Fisher Scientific Co., 2000 Park Ln., Pittsburg, PA 15275 or equivalent

²³Cat. No. 3246-01, J.T. Baker, Inc. or equivalent

²⁴Cat. No. M33-500, Fisher Scientific Co. or equivalent

Supplemental Assay Method for the Titration of Canine Parvovirus
in Cell Culture

2.2.6.6 Dissolve 1.03 g sodium phosphate, dibasic, anhydrous (Na_2HPO_4)²⁵ with 10 ml DW, heat to $60^\circ \pm 2^\circ\text{C}$ until dissolved; add to **Section 2.2.6.5** with constant mixing.

2.2.6.7 Dissolve 0.1 g calcium chloride, anhydrous (CaCl_2)²⁶ with 10 ml DW; add slowly to **Section 2.2.6.6** to avoid precipitation.

2.2.6.8 Q.S. to 1000 ml with DW, adjust pH to 7.0-7.3 with 2N HCl.

2.2.6.9 Sterilize through a 0.22- μm filter.

2.2.7 Glass slides, 8 chamber²⁷ (Lab-Tek® Slides)

2.2.8 Polystyrene tubes, 12 x 75 mm²⁸

2.2.9 Goat anti-cat IgG (H&L) fluorescein isothiocyanate labeled conjugate (Anti-cat Conjugate)²⁹

2.2.10 100% Acetone³⁰

2.2.11 Tuberculin syringe,³¹ 1 ml and needle, 20 ga x 1.5 in³²

2.2.12 Self-refilling repetitive syringe, 2 ml³³

2.2.13 Pipette, 10 ml³⁴

²⁵Cat. No. 3828-01, J.T. Baker, Inc. or equivalent

²⁶Cat. No. 4225-05, J.T. Baker, Inc. or equivalent

²⁷Cat. No. 177402, Nunc, Inc., 2000 N. Aurora Rd., Naperville, IL 60563 or equivalent

²⁸Falcon® 2058, Becton Dickinson Labware, 1 Becton Dr., Franklin Lakes, NJ 07417 or equivalent

²⁹Cat. No. 102-015-003, Jackson ImmunoResearch Laboratories, Inc., 872 W. Baltimore Pike, West Grove, PA 19390 or equivalent

³⁰Cat. No. A18-4, Fisher Scientific Co. or equivalent

³¹Cat. No. 309602, Becton Dickinson Labware or equivalent

³²Cat. No. 250107, Becton Dickinson Labware or equivalent

³³Wheaton®, Cat. No. 13-689-50C, Fisher Scientific Co. or equivalent

³⁴Falcon 7530, Becton Dickinson Labware or equivalent

Supplemental Assay Method for the Titration of Canine Parvovirus
in Cell Culture

3. Preparation for the test

3.1 Personnel qualifications/training

Personnel shall have experience in the preparation and maintenance of cell culture as well as in the propagation of animal viruses.

3.2 Preparation of equipment/instrumentation

3.2.1 On the day of inoculation, set a water bath at $36^{\circ} \pm 2^{\circ}\text{C}$.

3.2.2 On the day of the IFA test, prepare a humidity chamber in the aerobic incubator by filling a pan in the bottom with DW.

3.3 Preparation of reagents/control procedures

3.3.1 Preparation of CRFK Slides

3.3.1.1 Cells are prepared from healthy, confluent CRFK cells, that are maintained by passing every 3 to 4 days. On the day of test initiation, add 0.4 ml/chamber of $10^{5.2}$ to $10^{5.4}$ cells/ml diluted in Growth Media into all chambers of a Lab-Tek® slide. Prepare sufficient Lab-Tek® Slides to allow 25 wells for controls and 20 wells for each Test Serial. Incubate at $36^{\circ} \pm 2^{\circ}\text{C}$ in a CO_2 incubator. These become the CRFK Slides.

3.3.1.2 Use seeded CRFK Slides within 4 hr.

3.3.2 Preparation of the CPV Working Reference

3.3.2.1 On the day of inoculation, rapidly thaw a vial of CPV Reference in a $36^{\circ} \pm 2^{\circ}\text{C}$ water bath.

3.3.2.2 With the self-refilling repetitive syringe, dispense 1.8 ml of MEM into each of 7, 12 x 75-mm polystyrene tubes labeled 10^{-1} through 10^{-7} .

3.3.2.3 With a micropipettor, transfer 200 μl of the CPV Reference to the tube labeled 10^{-1} ; mix by vortexing.

**Supplemental Assay Method for the Titration of Canine Parvovirus
in Cell Culture**

3.3.2.4 Using a new pipette tip, transfer 200 μ l from the 10^{-1} labeled tube (**Section 3.3.2.3**) to the 10^{-2} tube; mix by vortexing.

3.3.2.5 Repeat **Section 3.3.2.4** for each subsequent dilution, transferring 200 μ l from the previous dilution to the next dilution tube until the tenfold dilution series is completed.

3.3.3 Preparation of Working FPV Antiserum

On the day of the IFA test, dilute FPV Antiserum in DPBS to the IFA working dilution on the CVB-L Reference and Reagent Sheet or as determined for that specific antiserum.

3.3.4 Preparation of Working Anti-cat Conjugate

On the day of the IFA test, dilute Anti-cat Conjugate in DPBS to the working dilution according to the manufacturer's recommendations.

3.4 Preparation of the sample

3.4.1 The initial test of a Test Serial will be with a single vial (a single sample from 1 vial). On the day of inoculation, rehydrate a vial of the Test Serial by transferring 1.0 ml for a 1-ml-dose vaccine, 0.5 ml for 1/2-ml-dose vaccines, etc., of the provided diluent into the vial containing the lyophilized Test Serial. Use a sterile 1.0-ml syringe and an 18-ga x 1.5-in needle; mix by vortexing. Incubate for 15 ± 5 min at room temperature (RT) ($23^{\circ} \pm 2^{\circ}\text{C}$).

3.4.2 Using the self-refilling repetitive syringe, dispense 1.8 ml MEM into each of 6, 12 x 75-mm polystyrene tubes labeled 10^{-1} through 10^{-6} .

Supplemental Assay Method for the Titration of Canine Parvovirus
in Cell Culture

3.4.3 With a micropipettor, transfer 200 μ l from the rehydrated Test Serial vial to the tube labeled 10^{-1} ; mix by vortexing.

3.4.4 Using a new pipette tip, transfer 200 μ l from the 10^{-1} labeled tube (**Section 3.4.3**) to the 10^{-2} tube; mix by vortexing.

3.4.5 Repeat **Section 3.4.4** for each subsequent dilution, transferring 200 μ l from the previous dilution to the next dilution tube, until the tenfold dilution series is completed.

4. Performance of the test

4.1 Inoculate 5 chambers/dilution of a CRFK Slide with 100 μ l/chamber from dilutions 10^{-6} through 10^{-3} of the Test Serial. Tip changes are not necessary between each dilution in a series if pipetting from the most dilute to the most concentrated within that series (e.g., 10^{-6} through 10^{-3}).

4.2 Inoculate 5 chambers/dilution of a CRFK Slide with 100 μ l/chamber, from dilutions 10^{-7} through 10^{-4} of the CPV Working Reference.

4.3 Five uninoculated chambers serve as a Negative Cell Control.

4.4 Incubate CRFK Slides in a $36^{\circ} \pm 2^{\circ}\text{C}$ CO_2 incubator for 120 ± 12 hr.

4.5 Following incubation, decant the media from the CRFK Slide and remove the plastic walls by twisting them away from the CRFK Slide, leaving the gasket attached to the slide.

4.6 Place the CRFK Slides in a slide rack; place the rack in a glass staining dish filled with DPBS. Incubate for 15 ± 5 min at RT.

4.7 Discard the DPBS; fix the CRFK Slides in 100% Acetone for 15 ± 5 min at RT. Remove and allow to air dry.

**Supplemental Assay Method for the Titration of Canine Parvovirus
in Cell Culture**

4.8 Pipette 75 ± 25 μ l of the Working FPV Antiserum into each chamber of the CRFK Slides. Incubate for 30 ± 5 min in the aerobic incubator at $36^{\circ} \pm 2^{\circ}\text{C}$.

4.9 Wash per **Section 4.6**.

4.10 Pipette 75 ± 25 μ l of the Working Anti-cat Conjugate into each chamber of the CRFK Slides. Incubate for 30 ± 5 min in the aerobic incubator at $36^{\circ} \pm 2^{\circ}\text{C}$.

4.11 Wash per **Section 4.6**.

4.12 Rinse the CRFK Slides with DW; allow to air dry.

4.13 Read at 100-200X magnification with a UV-light microscope; examine the cell monolayer for typical CPV apple-green nuclear fluorescence.

4.13.1 Chambers containing 1 or more cells with specific CPV fluorescence are considered to be positive.

4.13.2 Results are recorded as number of CPV positive chambers versus total number of chambers examined for each dilution of a Test Serial and the CPV Working Reference.

4.14 Calculate the 50% tissue culture infective dose (TCID_{50}) of the Test Serial and the CPV Working Reference using the method of Spearman-Kärber as modified by Finney.

Example:

10^{-3} dilution of Test Serial = 5/5 chambers IFA positive
 10^{-4} dilution of Test Serial = 5/5 chambers IFA positive
 10^{-5} dilution of Test Serial = 2/5 chambers IFA positive
 10^{-6} dilution of Test Serial = 0/5 chambers IFA positive

Spearman-Kärber calculation of total IFA positive wells (12), using 5 wells per dilution = 1.9 log

Log₁₀ of reciprocal of dose factor:

Titer of the Test Serial is $10^{5.9}$ TCID₅₀.

5.1.1 The calculated TCID₅₀ titer of the CPV Working Reference must fall within ± 2 standard deviations (SD) of its mean titer, as established from a minimum of 10 previously determined titers.

5.1.2 The lowest inoculated dilution of the CPV Working Reference must exhibit a 100% positive IFA reaction (5/5), and the highest (most dilute) must exhibit a negative IFA reaction (0/5).

5.1.3 The Negative Cell Control must not exhibit any cytopathic effect, specific CPV fluorescence, or cloudy media that would indicate contamination.

5.2 If the validity requirements are not met, then the assay is considered a **NO TEST** and may be retested without prejudice.

5.3 If the validity requirements are met and the titer of the Test Serial is greater than or equal to the titer contained in the Animal and Plant Health Inspection Service (APHIS) filed Outline of Production for the product under test, the Test Serial is considered **SATISFACTORY**.

Supplemental Assay Method for the Titration of Canine Parvovirus
in Cell Culture

5.4 If the validity requirements are met, but the titer of the Test Serial is less than the required minimum titer contained in the APHIS filed Outline of Production for the product under test, the Test Serial may be retested in accordance with the 9 CFR, Part 113.8.

6. Report of test results

Report results as the TCID₅₀ per dose of the Test Serial.

7. References

7.1 Code of Federal Regulations, Title 9, Part 113.317, U.S. Government Printing Office, 1999.

7.2 Cottral, GE, (Ed.), 1978. *Manual of standardized methods for veterinary microbiology*. Comstock Publishing Associates, Ithaca, NY. pg. 731.

7.3 Finney, DJ, 1978. *Statistical methods in biological assay*. Griffin, London. 3rd edition, pg. 508.

8. Summary of revisions

This document was rewritten to meet the current NVSL/CVB QA requirements, to clarify practices currently in use in the CVB-L, and to provide additional detail. No significant changes were made from the previous protocol.