

Center for Veterinary Biologics  
and  
National Veterinary Services Laboratories  
Testing Protocol

Supplemental Assay Method for Titrating the Fractions  
of Combination Avian Encephalomyelitis/Pox Vaccine

Date: December 16, 1998

Supersedes: June 1, 1984

Number: PYSAM0409.01

Standard Requirement: 9 CFR 113.325 and 113.326

Contact Person: Rebecca L. W. Hyde, (515) 239-8336

Approvals:

\_\_\_\_\_  
Rebecca L. W. Hyde, Head  
Poultry Biologics Virology Section  
Date: \_\_\_\_\_

\_\_\_\_\_  
P. Frank Ross, Acting Quality Assurance Manager  
Date: \_\_\_\_\_

\_\_\_\_\_  
\_/\_s/ Randall L. Levings\_\_\_\_\_  
Randall L. Levings, Director  
Center for Veterinary Biologics-Laboratory  
Date: \_12/16/98

United States Department of Agriculture  
Animal and Plant Health Inspection Service  
P. O. Box 844  
Ames, IA 50010

Mention of trademark or proprietary product does not constitute a guarantee or warranty of the product by USDA and does not imply its approval to the exclusion of other products that may be suitable.

Supplemental Assay Method for Titrating the Fractions of Combination Avian  
Encephalomyelitis/Pox Vaccine

Table of Contents

1. Introduction
2. Materials
  - 2.1 Equipment/instrumentation
  - 2.2 Reagents/supplies
3. Preparation for the test
  - 3.1 Personnel qualifications/training
  - 3.2 Preparation of equipment/instrumentation
  - 3.3 Preparation of reagents/control procedures
  - 3.4 Preparation of the sample
4. Performance of the test
  - 4.1 AE
  - 4.2 Pox
5. Interpretation of the test results
  - 5.1 Controls
  - 5.2 Calculating the titer
  - 5.3 Retests
  - 5.4 Evaluation of test results
6. Report of test results
7. References
8. Changes

Supplemental Assay Method for Titrating the Fractions of Combination Avian  
Encephalomyelitis/Pox Vaccine

## 1. Introduction

This Supplemental Assay Method (SAM) describes a procedure for titrating a vaccine containing both avian encephalomyelitis (AE) vaccine virus and avian pox vaccine virus. The vaccine is composed of separate preparations of each virus that are mixed together with a suitable stabilizer.

## 2. Materials

### 2.1 Equipment/instrumentation

Equivalent equipment or instrumentation may be substituted for any brand name listed below.

2.1.1 Centrifuge (Beckman J-6B, JS-4.2 rotor)

2.1.2 Humidified, rotating egg incubator (Midwest Incubators, Model No. 252)

2.1.3 Vortex mixer (Thermolyne Maxi Mix II, Model No. M37615)

2.1.4 Microliter pipette (Rainin Pipetman, P1000, or equivalent)

2.1.5 Cool-lite tester (Val-A)

2.1.6 Egg candling light on stand (Speed King)

2.1.7 Etcher electric engraver (Vibro-graver Acme Burgess, Inc.)

2.1.8 Glass 50-ml centrifuge tube, sterile (Kimax, Cat. No. 45166)

### 2.2 Reagents/supplies

Equivalent reagents or supplies may be substituted for any brand name listed below. All reagents and supplies must be sterile.

Supplemental Assay Method for Titrating the Fractions of Combination Avian  
Encephalomyelitis/Pox Vaccine

- 2.2.1 Cotton swabs/cotton balls
- 2.2.2 Serological pipettes (Falcon, Cat. No. 7530)
- 2.2.3 Pipette tips (Rainin Clean-Pak disposable microliter pipette tips RT-200)
- 2.2.4 Syringe, Glaspak 1 cc tuberculin, frosted tip, single use (Becton, Dickinson and Company)
- 2.2.5 Hypodermic needle, 26 ga x  $\frac{3}{8}$  in (Becton, Dickinson and Company, PrecisionGlide needle)
- 2.2.6 Hypodermic needle, 22 ga x  $1\frac{1}{2}$  in (Becton, Dickinson and Company, PrecisionGlide needle)
- 2.2.7 Glass test tubes, 16 x 125 with morton closures
- 2.2.8 Duco cement
- 2.2.9 Chloroform, reagent grade
- 2.2.10 Chick embryos from specific-pathogen-free (SPF) source
  - 1. Use 5- to 6-day-old embryos for the AE titration.
  - 2. Use 9- to 11-day-old embryos for the pox titration.
- 2.2.11 Solutions
  - 1. Tryptose Phosphate Broth (TPB)

TPB	29.5 g
q.s. with distilled or deionized water	1000.0 ml

Sterilize by autoclaving

Supplemental Assay Method for Titrating the Fractions of Combination Avian  
Encephalomyelitis/Pox Vaccine

**2. Penicillin/Streptomycin (pen/strep)**

penicillin g	15.775	g
streptomycin	100.0	g
q.s. with distilled or deionized water	1000.0	ml

Sterilize by filtration

**3. Normal Saline**

NaCl	8.5	g
q.s. with distilled or deionized water	1000.0	ml

Sterilize by autoclaving

**4. 70% alcohol**

ethyl alcohol	70	ml
q.s. with distilled or deionized water	30	ml

**5. Iodine, 2% in alcohol**

iodine	2	g
ethyl alcohol (70%)	100	ml

**2.2.12 Sterile distilled or deionized water**

**3. Preparation for the test**

**3.1 Personnel qualifications/training**

The executor must have experience or training in this protocol. This includes knowledge of aseptic biological laboratory techniques and preparation, proper handling and disposal of biological agents, reagents, tissue culture samples, and chemicals. The executor must also have knowledge of safe operating procedures and policies and Quality Assurance (QA) guidelines of the Center for Veterinary Biologics-Laboratory (CVB-L) or equivalent; and training in the operation of the necessary laboratory equipment listed in **part 2.1**.

Supplemental Assay Method for Titrating the Fractions of Combination Avian  
Encephalomyelitis/Pox Vaccine

### 3.2 Preparation of equipment/instrumentation

Operate all equipment/instrumentation according to manufacturers' instructions, and monitor in compliance with current corresponding CVB-L/National Veterinary Services Laboratories (NVSL) Standard Operating Procedures (SOPs) or equivalent.

### 3.3 Preparation of reagents/control procedures

Prepare reference viruses in the same manner as sample preparation.

### 3.4 Preparation of the sample

#### 3.4.1 AE

Rehydrate 500 doses of vaccine with 10.0 ml sterile purified water. Mix thoroughly. Transfer 1.0 ml of this vaccine to a 50-ml sterile glass centrifuge tube containing 9.0 ml sterile purified water. Mix thoroughly. Add 10.0 ml chloroform. Mix on a Vortex mixer for 3 separate 30-sec intervals. Centrifuge at 600 X *g* for 10 min. The aqueous phase is considered the  $10^0$  concentration of virus and contains 1 dose in 0.2 ml. Transfer 0.5 ml of the aqueous phase (upper layer) to a test tube containing 4.5 ml sterile purified water. Make further tenfold dilutions through  $10^{-5}$  using sterile purified water.

#### 3.4.2 Pox

Rehydrate 500 doses of vaccine with 10.0 ml sterile purified water. Mix thoroughly. Transfer 0.5 ml of this vaccine to a sterile test tube containing 4.5 ml of TPB. This is considered the  $10^0$  concentration and contains 1 dose in 0.2 ml. Make further tenfold dilutions, transferring 0.5 ml vaccine to 4.5 ml TPB, up through  $10^{-6}$ .

Supplemental Assay Method for Titrating the Fractions of Combination Avian  
Encephalomyelitis/Pox Vaccine

#### **4. Performance of test**

##### **4.1 AE**

Inoculate dilutions  $10^{-1}$  through  $10^{-5}$  into the yolk sac using 10 embryos per dilution. Inoculate 0.2 ml per embryo. Also have 20 uninoculated negative controls. Incubate the embryos and calculate the titer according to the criteria specified in the Code of Federal Regulations, Title 9 (9 CFR) 113.325(c)(2)(i).

##### **4.2 Pox**

Inoculate dilutions  $10^{-2}$  through  $10^{-6}$  onto the dropped chorioallantoic membrane (CAM) using at least 6 embryos per dilution. Inoculate 0.2 ml per embryo. Incubate the embryos and calculate the titer according to the criteria specified in 9 CFR 113.326.

#### **5. Interpretation of the test results**

##### **5.1 Controls**

Titrate a known positive reference virus with each group of titrations. The titer of the positive reference virus must be within the established range for the test results to be valid.

##### **5.2 Calculating the titer**

Determine the log 10 EID<sub>50</sub> titer using the method of Reed and Muench. This dilution and inoculation procedure allows for a direct readout on a per-dose basis. Round to 1 decimal.

##### **5.3 Retests**

Conduct retests as required by 9 CFR 113.8(b) and requirements of minimum release in firm's current Outline of Production, Part V.

Supplemental Assay Method for Titrating the Fractions of Combination Avian  
Encephalomyelitis/Pox Vaccine

#### 5.4 Evaluation of test results

5.4.1 The 9 CFR 113.8(b) defines the criteria for a satisfactory/unsatisfactory serial.

5.4.2 The firm's requirements of minimum release/stability titers for each vaccine are listed in the current Outline of Production, Part V, for the specific product code.

#### 6. Report of test results

Titers are reported out as EID<sub>50</sub> per bird dose.

#### 7. References

7.1 Reed, L.J., and H. Muench. 1938. A simple method of estimating 50% endpoints. Am. J. Hyg. 27:493-497.

#### 8. Changes

8.1 This document was rewritten to meet the current CVB-L QA SAM format. No significant changes were made from the previous protocol. This document supersedes the June 1, 1984, version.