

Center for Veterinary Biologics
and
National Veterinary Services Laboratories
Testing Protocol

Supplemental Assay Method for the Mikrozone
Electrophoretic Separation of Blood Fractions Used
in Veterinary Biologics

Date: April 24, 1997

Supersedes: TCSAM0509.01

Number: TCSAM0509.02

Contact Person: Jan DeJong, 515-239-8511
email: jdejong@aphis.usda.gov

Approvals:

P. Frank Ross, Head
Toxicology and Chemistry Section

Date:_____

Arthur J. Davis, Chief
Pathobiology Laboratory

Date:_____

Richard E. Hill, Quality Assurance Manager

Date:_____

_/_s/ R. L. Levings
Randall L. Levings, Director
Center for Veterinary Biologics-Laboratory

Date:_5/21/97_

United States Department of Agriculture
Animal and Plant Health Inspection Service
P. O. Box 844
Ames, IA 50010

Mention of trademark or proprietary product does not constitute a guarantee or warranty of the product by USDA and does not imply its approval to the exclusion of other products that may be suitable.

Supplemental Assay Method for the Mikrozone Electrophoretic Separation of
Blood Fractions Used in Veterinary Biologics

Table of Contents

1. Introduction
2. Materials
 - 2.1 Equipment/instrumentation
 - 2.2 Reagent/supplies
3. Preparation for the test
 - 3.1 Personnel qualifications/training
 - 3.2 Preparation of equipment/instrumentation
 - 3.3 Preparation of reagents/control procedures
 - 3.4 Preparation of sample
4. Performance of the test
5. Interpretation of the test results
6. Report of test results
7. References
8. Appendices
 - 8.1 Typical Densitometer Integration of a Reference Serum

**Supplemental Assay Method for the Mikrozone Electrophoretic Separation of
Blood Fractions Used in Veterinary Biologics**

1. Introduction

This laboratory protocol describes the procedure used to determine the per cent albumin or globulin present in a serum, antiserum or antitoxin. Combined with an assay for total protein, it can determine the protien content of the serum.

2. Materials--Brand names are provided for reference only. Equivalent equipment or materials may be used.

2.1 Equipment/instrumentation

2.1.1 Electrophoresis equipment (microzone cell) with power supply--(Brinkman Sartorius Microphor)

2.1.2 Microzone densitometer

2.1.3 Cellulose acetate strips--Size 57 x 145 mm(Sartorius cellulose acetate for electrophoresis Catalog No. 11207)

2.1.4 Microzone sample applicator--(Beckman Catalog No. 324399)

2.1.5 Computer--(with Densitometer integration calculation program)

2.1.6 Common laboratory apparatus and glassware--beakers, broad-tipped forceps (tweezer style), white paper towels or white filter paper, pen with non-diffusing ink, stainless steel pans with lids (8.5x5x2), ruler, class A volumetric flasks, parafilm, pipets or pipetor with tips

2.2 Reagents/supplies--all chemicals reagent grade unless specified. Store at room temperature unless otherwise specified. Chemicals and solutions good for 1 year unless otherwise specified.

2.2.1 Tris-barbital-sodium barbital buffer--(Gelman High Resolution Buffer Cat# 51104), pH 8.8. Undiluted buffer stable for 2 years.

2.2.2 Ponceau S fixative dye solution--(Beckman Cat# 324340) Store stock at 4 C. Stock stable for 2 years.

2.2.3 Acetic Acid--(Glacial--greater than 95%)

**Supplemental Assay Method for the Mikrozone Electrophoretic Separation of
Blood Fractions Used in Veterinary Biologics**

2.2.4 Bovine serum reference--normal bovine serum [NVSL media# 40032, bovine serum (donor)]. Store frozen (-20 C).

2.2.5 Water--use deionized or distilled water, or water of equivalent purity.

3. Preparation for the test

3.1 Personnel--No specific training is required. Individual should have working knowledge of laboratory equipment listed in **Section 2**.

3.2 Equipment/instrumentation

3.2.1 Prepare electrophoresis cell

3.2.1.1 Place cell on horizontal surface and fill the buffer compartments to the fill line with tris-barbital-sodium barbital pH 8.8 buffer solution (**3.3.1**).

3.2.1.2 With laboratory toweling, wipe any buffer from the edges and center section of the cell.

3.2.2 Prepare cellulose acetate strips

Critical control point: Touch the cellulose acetate strip as little as possible, use tweezers.

3.2.2.1 Using a pen with non-diffusing ink, mark each strip with a reference line 2 inches from the left end of the strip. Write all sample numbers, dates and other information on the left end of the strip.

3.2.2.2 Gently float the strip on the surface of buffer solution which has been poured into a flat stainless steel pan. When the strip is wetted, use tweezers to gently submerge the strip in the buffer.

Critical control point: Do not trap air bubbles between the strip and the buffer solution. Be careful to not submerge the strip before it is thoroughly wetted.

**Supplemental Assay Method for the Mikrozone Electrophoretic Separation of
Blood Fractions Used in Veterinary Biologics**

3.2.2.3 After approximately 5 minutes, carefully remove strip from the buffer and remove excess buffer by gently blotting it between two sheets of white toweling or two filter paper sheets.

Critical control point: Blot only until excess moisture is removed--Do not over dry (dry spots will appear on the strip).

3.2.2.4 Carefully stretch the blotted cellulose acetate strip across the bracket and place in the electrophoresis cell with the ends of the strip hanging directly in the buffer.

3.3 Reagents/controls

3.3.1 Tris-barbital-sodium barbital buffer--dilute contents of one packet with water following label directions. Mix until completely dissolved. Buffer stable for 6 to 8 weeks if stored at 4 C, warm to room temperature before using.

3.3.2 Ponceau S working solution--dilute contents of one container of stock dye with water following label directions.

3.3.3 Acetic acid, 5%--carefully add 50 ml glacial acetic acid to 600 ml water in 1000 ml volumetric flask. Mix and fill to volume with water.

3.3.4 Serum reference--thaw serum, warm to room temperature and mix gently by inverting vial 6--8 times.

3.4 Sample

3.4.1 Serum products (blood origin products) and serum samples are usually ready for electrophoresis. Serum fractions produced by ammonium sulfate or sodium sulfate precipitation must be dialyzed to remove high salt content. See current version of **TCSOP0001** for receipt of sample.

3.4.1.1 Using a pipet or pipetor with tip, place a drop of sample (gently mixed) on a sheet of parafilm with a sample number marked to identify each drop. (Include reference serum with each run.)

**Supplemental Assay Method for the Mikrozone Electrophoretic Separation of
Blood Fractions Used in Veterinary Biologics**

3.4.1.2 Use a microzone sample applicator to apply the samples to the strip, with the application points near the reference line (approximately 1/2 inch to the right of the line). Be careful that only the tip of the applicator touches the sample droplet and that the applicator tip is rinsed and blotted dry between sample applications. Run duplicates of each sample and the reference serum. (Run duplicates of the reference on each strip.)

4. Test Procedure

4.1 Electrophoresis run

4.1.1 Place cover on the electrophoresis cell and, attach electrical connections to the power supply. Run current from "left" to "right" end of strip.

4.1.2 Set current to constant voltage and switch power on. Adjust voltage to 250 volts.

4.1.3 Run migration for 25 minutes.

4.1.4 Turn off power and disconnect cell from the power supply. The same batch of buffer may be used for more than one run but the current must be reversed with each additional run (changing the sides of the anode and the cathode on the cell). The strips must also be reversed, so migration is in the same direction on the strip.

4.2 Staining

4.2.1 Carefully remove the strip from the electrophoresis cell. Float strip on staining solution, as with buffer (3.2.2.2), and then submerge. Allow to stain for approximately 5 minutes.

Critical control point: Do not place one strip over another in the staining solution, stain strips individually.

4.2.2 Remove strip from the stain and rinse in 3 rinses of 5% acetic acid--approximately 5 minutes in each rinse. If the background is not white, repeat the rinses.

**Supplemental Assay Method for the Mikrozone Electrophoretic Separation of
Blood Fractions Used in Veterinary Biologics**

4.2.3 Remove the strip from the last rinse and blot with clean white toweling. Check to determine if the protein fractions appear as distinct red bands. If bands appear in proper places and stain well, place the strip flat between sheets of toweling and press with a heavy, broad-surfaced object until dry.

4.2.4 If reference sample does not migrate or stain properly, repeat run.

5. Calculations/evaluation

5.1 Integrate patterns with the densitometer. Use the ink mark on the strip as an end point reference. Determine the albumin and globulin percentages using the area of the peaks in the integrated pattern.

5.2 Determine the albumin (A) or globulin (G) fraction (%) by dividing the A (or G) integrator area by the total integrator area. To determine the A (or G) concentration (mg/ml), multiply that fraction by the total protein value of the sample (determine total protein concentration according to current version of **TCSAM0504**). A computer program can be used for these calculations (In the TC Section, a computer program called **Electrophoresis** resident on the computer in W6 is used for these calculations). The A/G ratio is normally calculated for each sample. A typical densitometer integration of a reference serum is shown in **Appendix 8.1**.

6. Report of test results

6.1 Test results are reported following current version **TCSOP0001**

7. References

7.1 Reference is made to the long time use of this procedure by the Toxicology/Chemistry Section.

7.2 Version .01 was generated as a revision of SAM 509, December 15, 1970, to meet format requirements. Version .02 represents format changes resulting from the creation of the Center for Veterinary Biologics. No technical changes from the previous version are represented.

**Supplemental Assay Method for the Mikrozone Electrophoretic Separation of
Blood Fractions Used in Veterinary Biologics**

7.3 For a hard copy of the appendix, please contact Jan
DeJong 515-239-8511.

8. Appendices

8.1 Typical Densitometer Integration of a Reference Serum