

**Center for Veterinary Biologics
and
National Veterinary Services Laboratories
Testing Protocol**

**Supplemental Assay Method for the Determination of
Protein and Phenol in PPD (Purified Protein Derivative
Produced From Cultures of *Mycobacterium bovis* Strain
AN-5) Bovis Tuberculin**

Date: December 10, 1998
Supersedes: TCSAM0513.03, April 1, 1998
Number: TCSAM0513.04
Standard Requirement: 9 CFR 113.409
Contact Person: Gary B. Clark, (515) 239-8514
E-mail: Gary.B.Clark@usda.gov
Approvals:

P. Frank Ross, Head
Toxicology and Chemistry Section
Date: _____

Arthur J. Davis, Chief
Pathobiology Laboratory
Date: _____

P. Frank Ross, Acting Quality Assurance Manager
Date: _____

_/_s/ Randall L. Levings_____
Randall L. Levings, Director
Center for Veterinary Biologics-Laboratory
Date: _12/11/98

_/_s/ J. E. Pearson_____
James E. Pearson, Director
National Veterinary Services Laboratories
Date: _12/11/98

United States Department of Agriculture
Animal and Plant Health Inspection Service
P. O. Box 844
Ames, IA 50010

Mention of trademark or proprietary product does not constitute a guarantee or warranty of the product by USDA and does not imply its approval to the exclusion of other products that may be suitable.

Supplemental Assay Method for the Determination of Protein and Phenol in PPD
(Purified Protein Derivative Produced From Cultures of *Mycobacterium bovis* Strain
AN-5) Bovis Tuberculin

Table of Contents

1. Introduction
2. Materials
 - 2.1 Equipment
 - 2.2 Reagents/supplies
 - 2.2.1 Protein Test
 - 2.2.2 Phenol Test
3. Preparation for the test
 - 3.1 Personnel qualifications/training
 - 3.2 Preparation of equipment/instrumentation
 - 3.3 Preparation of reagents/control procedures
 - 3.3.1 Protein test
 - 3.3.2 Phenol test
 - 3.4 Preparation of the sample
 - 3.4.1 Receipt
 - 3.4.2 Preparation
4. Performance of the test
 - 4.1 Protein
 - 4.2 Phenol
5. Interpretation of the test results
 - 5.1 Protein
 - 5.2 Phenol
 - 5.3 Controls
6. Report of test results
7. References
8. Appendices
 - 8.1 TC Log--Tuberculin (PPD) Log
9. Quick reference

Supplemental Assay Method for the Determination of Protein and Phenol in PPD
(Purified Protein Derivative Produced From Cultures of *Mycobacterium bovis* Strain
AN-5) Bovis Tuberculin

1. Introduction

The Code of Federal Regulations, Title 9 (9 CFR) (Animals and Animal Products) states that the Animal and Plant Health Inspection Service (APHIS) is responsible for administering the Virus-Serum-Toxin Act. It requires licensed tuberculin be tested and approved by the Center for Veterinary Biologics (CVB) before it can be marketed. Protein concentration is determined by classical Kjeldhal digestion, distillation, and titration of the ammonia. Phenol is determined by end-point titration with bromate/bromide. Satisfactory product must contain 1.0 mg/ml \pm 0.1 mg/ml protein. Phenol content must be 0.50% \pm 0.04%.

2. Materials

2.1 Equipment

2.1.1 Balance, top loading, capable of measuring 0.01 g

2.1.2 Digestion unit, Buchi, B-426, with digestion tubes

2.1.3 Distillation unit, Buchi, B-316

2.1.4 Volumetric pipets, Class A, meets ASTM Standard E969-83

2.1.5 Volumetric flasks, Class A, with barrel head glass stopper, meets ASTM E288 requirements

2.1.6 125-ml erlenmeyer flasks

2.1.7 10-ml buret with PTFE stopcock, precision bore, calibrated to ASTM E-694 accuracy requirements

2.1.8 50-ml buret with PTFE stopcock, precision bore, calibrated to ASTM E-694 requirements

2.1.9 Graduated cylinders, 50, 100, 250, 500, and 1,000 ml, (PYREX), meets ASTM D86, D216, D447 requirements

2.1.10 Glass-stoppered erlenmeyer flasks, 250 ml

Supplemental Assay Method for the Determination of Protein and Phenol in PPD
(Purified Protein Derivative Produced From Cultures of *Mycobacterium bovis* Strain
AN-5) Bovis Tuberculin

2.1.11 Heating/stirring plate with stirring bars

2.1.12 Fast filter paper

2.2 Reagents/Supplies: All chemicals, reagent grade. Use distilled or demineralized water or water of equivalent purity.

2.2.1 Protein

1. Sulfuric acid (H_2SO_4)--Purity: Minimum 95.0%, Maximum 98.0%
2. Mercury Tablets, Brinkmann Instruments, Catalog No. 015-00-646-3
3. Sodium hydroxide (NaOH)--Purity: 98.5%
4. Boric acid (H_3BO_3)--Purity: 99.9%
5. Methyl red--Purity: 98.0%
6. Hydrochloric acid (HCl)--Assay: 36.5%-38.0%
7. Sodium carbonate (Na_2CO_3)--Purity: 99.9%
8. Bromo phenol blue--Purity: 98.0%
9. National Veterinary Services Laboratories (NVSL) Control--Pool of PPD Tuberculin products with established protein and phenol values
10. Protein, Bovuminar crystallized powder--Intergen Company, 2 Manhattanville Road, Purchase, NY 10577, Catalog No. 3000-70, Purity: 98.9%

2.2.2 Phenol (some reagents same as for protein)

1. Methyl orange--Purity: 98.0%
2. Silicotungstic acid
($\text{H}_4[\text{Si}(\text{W}_3\text{O}_{10})_4] \cdot 26\text{H}_2\text{O}$)--Purity: 99.0%

Supplemental Assay Method for the Determination of Protein and Phenol in PPD
(Purified Protein Derivative Produced From Cultures of *Mycobacterium bovis* Strain
AN-5) Bovis Tuberculin

3. Arsenic trioxide (As_2O_3)--Purity: 99.9%
4. Sodium bicarbonate (NaHCO_3)--Purity: 99.9%
5. Potassium bromate (KBrO_3)--Purity: 98.5%
6. Potassium bromide (KBr)--Purity: 99.0%
7. Phenol ($\text{C}_6\text{H}_5\text{OH}$)--Purity: $\geq 99.0\%$

3. Preparation for the Test

3.1 Training of technical personnel

No special test-related training is needed for this testing. Analysts performing this procedure should first conduct 2 trial runs using controls and standards and obtain results within acceptable limits.

3.2 Preparation of equipment/instrumentation

Become familiar with Buchi instruction regarding operation. Turn on water that aspirates fume from suction tube of the digestion unit and keep the water cool in the condenser of the distillation unit. Adjust water flow in the distillation unit to approximately 1 L per min. Turn on the distillation unit. Set time preselector to "2" (2 min) and stopcock for aspiration to "Off." Make sure that Buchi bottles of NaOH and water are adequate.

3.3 Preparation of reagents

3.3.1 Protein Test (all reagents stable for at least 6 mo unless specified)

1. Cut Hg tablets into half.

Caution: Because tablet contains mercury, handle in fume hood and wear gloves, protective glasses, and mask.

Supplemental Assay Method for the Determination of Protein and Phenol in PPD
(Purified Protein Derivative Produced From Cultures of *Mycobacterium bovis* Strain
AN-5) Bovis Tuberculin

2. 32% NaOH, dissolve 640 g \pm 1 g NaOH in 1.4 L H₂O in 2-L volumetric flask on the magnetic stirrer. Cool to room temperature. Dilute to volume with H₂O. Repeat above until Buchi 10-L bottle is full. Store at 4°C.

Caution: NaOH is caustic--Avoid contact with skin.

3. Saturated H₃BO₃, add 15 g to 100 ml H₂O. Stir, with heat, until all H₃BO₃ dissolves. Some H₃BO₃ recrystallizes when cool. Store at room temperature (RT).

4. 0.5% methyl red, dissolve 0.5 g in 100 ml ethanol. Store at RT. Stable for 2 mo.

5. Standardized 0.01 N HCl-0.02 N HCl, 1.7 ml HCl/L H₂O. Titrate exactly 0.0100 g dried sodium carbonate dissolved in 25 ml H₂O. Indicator: 3 drops 0.1% bromo phenol blue; the color of endpoint is green not bluish green, nor yellowish green. Store at RT.

Calculation:

$$\underline{N} \text{ HCl} = [(g \text{ Na}_2\text{CO}_3) \times (1000)] / [(Vol \text{ HCl}) \times (52.994)].$$

Caution: Concentrated HCl is corrosive--Handle in fume hood. Avoid contact with skin.

6. Protein standard, weigh approximately 2 g Bovuminar crystallized powder and transfer to 2-L flask. Dissolve and dilute to 2 L with H₂O. Aliquot into 15-ml portions in 30-ml serum vials. Seal under nitrogen. Store at 4°C.

3.3.2 Phenol Test (all reagents stable for at least 6 mo unless specified)

1. 20% HCl, slowly add 200 ml HCl to 600 ml H₂O, dilute to 1 L. Store at RT.

Supplemental Assay Method for the Determination of Protein and Phenol in PPD
(Purified Protein Derivative Produced From Cultures of *Mycobacterium bovis* Strain
AN-5) Bovis Tuberculin

2. 0.1% methyl orange, add 0.1 g methyl orange to 100 ml H₂O. Filter if necessary. Make fresh every 2 mo. Store at RT.

3. Silicotungstic acid solution (SAS), dissolve 60 g H₄[Si(W₃O₁₀)₄]*26H₂O in 400 ml H₂O in 500-ml volumetric flask. Add 50 ml H₂SO₄. When cool, dilute to volume with H₂O. Store at 4°C.

4. Clarifying solution (CS), add 50 ml SAS and 125 ml 20% HCl to 325 ml H₂O. Prepare fresh prior to each test.

5. "Acid solution" for As₂O₃ standard solution, add 110 ml HCl and 2.5 ml methyl orange to 100 ml H₂O. Store at RT.

6. 0.0500 N As₂O₃, dissolve 2.4730 g dried As₂O₃ in 25 ml hot 1N NaOH in 1-L volumetric flask. Neutralize with 25 ml 1N H₂SO₄. Cool and dilute to vol with H₂O. Store at RT.

Caution: As₂O₃ is extremely toxic--Avoid contact, handle in fume hood using gloves, mask, and goggles. Consult Material Safety Data Sheet for specific handling instructions.

7. Phenol standard, dissolve 0.50 g phenol in 50 ml H₂O and dilute to volume. Store at RT.

Critical Control Point: The final diluted volume of the test fluid must be adjusted as described in 4.3.2.8.

8. Test fluid (TF), dissolve 0.30 g NaHCO₃, 1.67 g KBrO₃, and 15.00 g KBr in H₂O and qs to 1 L with H₂O. Store at RT. The TF volume must be adjusted by adding corrected volume of H₂O to TF. It must take a volume of 21.3 ml to titrate 25 ml 0.050 N As₂O₃ in 10 ml "Acid Solution." A first time titration will require less than 21.3 ml TF. Adjust as described in the following example:

Supplemental Assay Method for the Determination of Protein and Phenol in PPD
(Purified Protein Derivative Produced From Cultures of *Mycobacterium bovis* Strain
AN-5) Bovis Tuberculin

Example: Assume the first time titration volume
is 20.5 ml.

$$(1,000 \text{ ml of TF}) - (20.5 \text{ ml}) = 979.5 \text{ ml}$$

$$\frac{(979.5)(\text{desired vol})}{(\text{actual vol})(20.5)} \text{ or } \frac{(979.5)(21.3)}{(20.5)} = 1,017.2 \text{ ml}$$

For corrected vol of H₂O: 1017.2 - 979.4 =
37.8 ml to be added to TF.

Note: TF in buret has to be put back into flask.

3.4 Preparation of the sample

3.4.1 Receipt--Reference current version of TCSOP0001.

3.4.2 Preparation--PPD tuberculin products are stored
at 4°C in the walk-in refrigerator prior to testing.
Before testing, allow sample vials and reagents to warm
to room temperature.

4. Performance of the test use Tuberculin (PPD) Log Sheet, Appendix 8.1

4.1 Protein--Analyze the control pool and protein standard
each time testing is performed. Analyze each in triplicate.

4.1.1 Place 5.0 ml PPD, ½ Hg tablet, and
3.0 ml ± 0.1 ml H₂SO₄ into a Buchi digestion tube.
Same for the standard and control.

Caution: HgO is poisonous--Handle in fume hood, using gloves,
mask, and goggles. Consult Material Safety Data Sheet for
specific handling instructions.

Caution: Concentrated H₂SO₄ is corrosive--Handle in fume
hood. Avoid contact with skin.

4.1.2 On digestion rack heater, digest for 30 min after
all H₂O is distilled and acid comes to true boil.

**Supplemental Assay Method for the Determination of Protein and Phenol in PPD
(Purified Protein Derivative Produced From Cultures of *Mycobacterium bovis* Strain
AN-5) Bovis Tuberculin**

4.1.3 Cool, add 3 ml H₂O and cool again. Add 15 ml NaOH solution to sample and attach to distillation apparatus.

4.1.4 Distill about 25 ml into 50-ml erlenmeyer flask containing 5 ml H₃BO₃ and 3 drops indicator.

4.1.5 Titrate to endpoint color change of yellow to deep rose (pH 5.0) with HCl. Record result on log sheet.

4.2 Phenol (Analyze the control pool and phenol standard each time testing is performed. Analyze each in triplicate.)

4.2.1 Add 5 ml PPD and 100 ml CS to 250-ml glass stoppered flask. Shake 2 min. Filter through fast filter paper.

4.2.2 Transfer 50 ml filtrate to another flask. Add 1 drop methyl orange, stopper, and shake a few sec. Observe the color, when red, go to **4.2.3**.

4.2.3 Titrate with 2 ml test fluid (TF), stopper, and shake a few sec and observe the color. When red, repeat **4.2.3**. When colorless, go to **4.2.4**.

4.2.4 Shake 30 sec. Add 1 drop indicator, stopper, and shake a few sec and observe the color. When it does not turn to colorless within 10 sec, titrate with 1 ml TF, stopper, and repeat **4.2.4**. When colorless, go to **4.2.5**.

4.2.5 Shake 1 min. Add 1 drop indicator, stopper, and shake a few sec. Observe the color. When red stays longer than 10 sec, titrate with 0.50-ml TF, stopper, and repeat **4.2.5**. When colorless, record total vol of TF as the endpoint of titration and use for calculation of percent phenol.

5. Interpretation of the test

5.1 Protein (Report average of triplicates.)

mg Protein/ml = (ml HCl)(N HCl)(1.4007)(6.25)/(5 ml PPD)
Satisfactory Protein Content: 1.0 mg/ml ± 0.1 mg/ml

5.2 Phenol (Report average of triplicates.)

Percent phenol = (vol of test fluid)(0.04)-(0.04)
Satisfactory Phenol Content: 0.50% ± 0.04%

Supplemental Assay Method for the Determination of Protein and Phenol in PPD
(Purified Protein Derivative Produced From Cultures of *Mycobacterium bovis* Strain
AN-5) Bovis Tuberculin

5.3 Controls

Results for controls and standards must be within acceptable limits; otherwise repeat testing.

6. Reporting of test results

Validate and report results according to the current version of TCSOP0001.

7. References

7.1 Code of Federal Regulations, Title 9, Section 113.409
Revised January 1992, page 629.

7.2 Official Methods of Analysis of AOAC International,
Arlington, Virginia, 16th Edition, Pat Cuniff, Editor
(1995), Volume I, Chapter 12, page 7.

8. Summary of Changes

Version .01 was written to meet NVSL/CVB Quality Assurance requirements, to clarify practices in use in the NVSL/CVB-L, and to provide additional detail. No significant changes were made from the previous protocol.

Version .02 was written to clarify practices in use in the NVSL/CVB-L and to provide additional detail.

Version .03 was written to clarify practices in use in the NVSL/CVB-L and to provide additional detail. The following are the significant changes made from the superseded protocol:

1. Change in the digestion apparatus

Version .04 was written to clarify practices in use in the NVSL/CVB-L and to provide additional detail.

Supplemental Assay Method for the Determination of Protein and Phenol in PPD
(Purified Protein Derivative Produced From Cultures of *Mycobacterium bovis* Strain
AN-5) Bovis Tuberculin

Appendix 8.1 TC Log Sheet

Tuberculin (PPD) Log

Cases #: _____	Date: _____
Sample Code: _____	Analyst: _____
Serial #: _____	Reviewed by: _____
Comment: _____	

Protein Micro-Kjedhal (5ml sample)		[Vol. HCl] [N HCl] [1.4008] [10] [6.25] (Vol. sample)	
Sample	N	mg Protein/ml	Average of 3
1 _____	a		_____ mg/ml
	b		
	c		
2 _____	a		_____ mg/ml
	b		
	c		
Std.: _____	a		_____ mg/ml
	b		
	c		
Control: _____	a		_____ mg/ml
	b		
	c		

Phenol Titration: (5ml sample)		[Vol. Test Fluid] [.04] - [.04]	
Sample	Vol. Test Fluid	% Phenol	Average of 3
1 _____	a		_____ %
	b		
	c		
2 _____	a		_____ %
	b		
	c		
Std.: _____	a		_____ %
	b		
	c		
Control: _____	a		_____ %
	b		
	c		

Supplemental Assay Method for the Determination of Protein and Phenol in PPD
(Purified Protein Derivative Produced From Cultures of *Mycobacterium bovis* Strain
AN-5) Bovis Tuberculin

9. Quick reference

- ____ Product code _____
- ____ Submission paperwork correct
- ____ Received correct specimens and adequate amount
- ____ Protein standard, date of origin _____
- ____ Standard acid
- ____ Test fluid, date of origin _____
- ____ Control PPD, date of origin _____
- ____ Triplicate determinations conducted
- ____ Control phenol result within limits
- ____ Control protein result within limits
- ____ Protein standard result within limits
- ____ Phenol standard result within limits
- ____ Results entered into computer, date _____
- ____ Results validated, date _____

Analyst _____ Date _____