## **Barcodes**

If barcodes are available from the lab, here is how to use them:

Barcodes are printed in sets of 8 individual labels (See figure 1 below) that should be used for one animal only. If you use barcodes you MUST attach them as follows:

- 1 label on the sample tube (As shown in figure 2)
- 1 label on the USDA BSE Surveillance Data Collection form (as shown in Figure 3)
- 1 label on the plastic bag containing ID devices (tattoos, brands etc) removed from the animal that was sampled

If appropriate they should also be used as follows:

- 1 label on the USDA BSE Surveillance Submission form in the field marked "BSE Sample ID"instead of writing in 1, 2, 3 etc. as noted in the instructions. (Because of size, you may not be able to get 20 bar code labels on the USDA BSE Surveillance Submission form – use the USDA BSE Surveillance Submission Continuation form for additional animals)
- 1 each on the head and carcass (attach to a bangle tag or other labeling device)
- 2 for any other applications as needed NOTE: Any labels that are not used should be destroyed



Figure 1



Figure 2



Figure 3

August 24, 2004 Barcodes: Page H1