

**Center for Veterinary Biologics  
and  
National Veterinary Services Laboratories  
Testing Protocol**

**Supplemental Assay Method for the Titration of  
Newcastle Disease Vaccine, Infectious Bronchitis  
Vaccine, and Combination Newcastle Disease-Infectious  
Bronchitis Vaccine in Chicken Embryos**

Date: October 29, 1998

Supersedes: September 1, 1985

Number: PYSAM0411.01

Standard Requirement: 9 CFR 113.327 and 113.329

Contact Person: Rebecca L. W. Hyde, (515) 239-8336

Approvals:

\_\_\_\_\_  
Rebecca L. W. Hyde, Head  
Poultry Biologics Virology Section

Date: \_\_\_\_\_

\_\_\_\_\_  
Richard E. Hill, Quality Assurance Manager

Date: \_\_\_\_\_

\_\_\_\_\_  
Randall L. Levings, Director  
Center for Veterinary Biologics-Laboratory

Date: \_\_\_\_\_

United States Department of Agriculture  
Animal and Plant Health Inspection Service  
P. O. Box 844  
Ames, IA 50010

Mention of trademark or proprietary product does not constitute a guarantee or warranty of the product by USDA and does not imply its approval to the exclusion of other products that may be suitable.

Supplemental Assay Method for the Titration of Newcastle Disease Vaccine,  
Infectious Bronchitis Vaccine, and Combination Newcastle Disease-Infectious  
Bronchitis Vaccine in Chicken Embryos

Table of Contents

1. Introduction
2. Materials
  - 2.1 Equipment/instrumentation
  - 2.2 Reagents/supplies
3. Preparation for the test
  - 3.1 Personnel qualifications/training
  - 3.2 Preparation of equipment/instrumentation
  - 3.3 Preparation of reagents/control procedures
  - 3.4 Preparation of the sample
4. Performance of the test
  - 4.1 Egg inoculation
  - 4.2 Incubation
5. Interpretation of the test results
  - 5.1 Controls
  - 5.2 NDV
  - 5.3 IBV
  - 5.4 Calculations
  - 5.5 Retests
  - 5.6 Evaluation of test results
6. Report of test results
7. References

**Supplemental Assay Method for the Titration of Newcastle Disease Vaccine,  
Infectious Bronchitis Vaccine, and Combination Newcastle Disease-Infectious  
Bronchitis Vaccine in Chicken Embryos**

## **1. Introduction**

This Supplement Assay Method (SAM) describes a procedure for titrating Newcastle disease (ND) vaccine, Infectious Bronchitis (IB) vaccine, and combination ND-IB vaccine. The vaccines are reconstituted and inoculated into embryonated chicken eggs in 10-fold dilutions such that the 50% endpoint of infectivity ( $EID_{50}$ ) can be calculated directly on a per field dose basis by the Reed-Muench method.

## **2. Materials**

### **2.1 Equipment/instrumentation**

Equivalent equipment or instrumentation may be substituted for any brand name listed below.

- 2.1.1 Centrifuge (Beckman J-6B, JS-4.2 rotor)
- 2.1.2 Humidified, rotating egg incubator (Midwest Incubators, Model 252)
- 2.1.3 Vortex mixer (Thermolyne Maxi Mix II Model No. M37615)
- 2.1.4 Pipette (Rainin Pipetman, P1000)
- 2.1.5 Cool-lite tester (Val-A)
- 2.1.6 Egg candling light on stand (Speed King)
- 2.1.7 Etcher electric engraver (Vibro-graver Acme Burgess, Inc.)

### **2.2 Reagents/supplies**

Equivalent reagents or supplies may be substituted for any brand name listed below. All reagents and supplies must be sterile.

- 2.2.1 Cotton swabs/cotton balls
- 2.2.2 Serological pipets (Falcon, Cat. No. 7530)
- 2.2.3 Specific-pathogen-free (SPF) chick embryos, 9- to 11-day-old

Supplemental Assay Method for the Titration of Newcastle Disease Vaccine,  
Infectious Bronchitis Vaccine, and Combination Newcastle Disease-Infectious  
Bronchitis Vaccine in Chicken Embryos

2.2.4 Disposable sterile filters, 0.45 um (Millex-HA Millipore)

2.2.5 Pipette tips (Rainin Clean-Pak disposable microliter pipette tips RT-200)

2.2.6 Syringe, glaspak, 1-cc tuberculin, frosted tip, single use (Becton, Dickinson and Company)

2.2.7 Hypodermic needle, 18 gaugex1½ in (Becton, Dickinson and Company, PrecisionGlide needle)

2.2.8 Hypodermic needle, 25 gaugex5⁄8 in (Becton, Dickinson and Company, PrecisionGlide needle)

2.2.9 Glass test tubes, 16x125 with morten closures

2.2.10 Glass test tubes, 13x100 with morten closures

2.2.11 Duco cement

2.2.12 Pipette tips (Rainin 0-100, 0-200, 100-1000 or equivalent)

2.2.13 Solutions

1. Tryptose phosphate broth (TPB)

TPB	29.5	g
q.s. with distilled or deionized water	1000.0	ml
Sterilize by autoclaving		

2. Penicillin/Streptomycin (pen/strep)

penicillin g	15.775	g
streptomycin	100.0	g
q.s. with distilled or deionized water	1000.0	ml

**Supplemental Assay Method for the Titration of Newcastle Disease Vaccine,  
Infectious Bronchitis Vaccine, and Combination Newcastle Disease-Infectious  
Bronchitis Vaccine in Chicken Embryos**

**3. Normal saline (0.85%)**

NaCl	8.5	g
q.s. with distilled or deionized water	1000.0	ml
Sterilize by autoclaving		

**4. 70% alcohol**

ethyl alcohol	70.0	ml
q.s. with distilled or deionized water	30.0	ml

**5. Iodine 2% in alcohol**

iodine	2.0	g
ethyl alcohol (70%)	100.0	ml

**3. Preparation for the test**

**3.1 Personnel qualifications/training**

The executor must have experience or training in this protocol. This includes knowledge of aseptic biological laboratory techniques and preparation, proper handling, and disposal of biological agents, reagents, tissue culture samples, and chemicals. The executor must also have knowledge of safe operating procedures and policies and Quality Assurance (QA) guidelines of the Center for Veterinary Biologics-Laboratory (CVB-L) or equivalent; and training in the operation of the necessary laboratory equipment listed in part 2.1.

**3.2 Preparation of equipment/instrumentation**

Operate all equipment/instrumentation according to manufacturer's instructions and monitor in compliance with current corresponding CVB-L/National Veterinary Services Laboratories Standard Operating Procedures (SOPs), or equivalent.

Supplemental Assay Method for the Titration of Newcastle Disease Vaccine,  
Infectious Bronchitis Vaccine, and Combination Newcastle Disease-Infectious  
Bronchitis Vaccine in Chicken Embryos

### 3.3 Preparation of reagents/control procedures

Prepare reference viruses in the same manner as sample preparation.

### 3.4 Preparation of the sample

#### 3.4.1 ND vaccine

1. Reconstitute ND vaccine, frozen or lyophilized, in 4°C tryptose phosphate broth (TPB) to a total vol as listed in **Table 1**. Vigorously mix for a minimum of 30 sec, and hold in an ice bath for 15 min to allow for virus disaggregation. Mix again and further dilute the reconstituted vaccine in TPB to obtain the  $10^0$  concentration (1 dose per 0.1 ml) as shown in **Table 1**.

2. Make 10-fold dilutions of the  $10^0$  virus concentration,  $10^{-1}$  through  $10^{-8}$ , by serially mixing 0.5 ml of the virus with 4.5 ml of TPB plus antibiotics (500-u penicillin and 2.0-mg streptomycin per ml TPB). Use a separate pipette for each transfer. Mix each dilution well before proceeding to the next dilution, and keep all dilutions on ice.

#### 3.4.2 IB vaccine

Prepare the same as for the ND vaccine (See **3.4.1**), making 10-fold dilutions through  $10^{-6}$ .

#### 3.4.3 Combination ND-IB vaccine

Prepare the same as for the ND vaccine. For titration of the Newcastle disease virus (NDV) fraction, the dilutions are inoculated without further treatment. For titration of the infectious bronchitis virus (IBV) fraction, neutralize the NDV by adding 0.8 ml of each dilution,  $10^{-2}$  through  $10^{-6}$ , to an equal vol (0.8 ml) of anti-NDV serum. Mix well, and incubate in an ice bath for 30 min prior to inoculation.

Supplemental Assay Method for the Titration of Newcastle Disease Vaccine,  
Infectious Bronchitis Vaccine, and Combination Newcastle Disease-Infectious  
Bronchitis Vaccine in Chicken Embryos

**Table 1**

Number of doses in 1 vial	Reconstituted Volume of Vaccine	Additional dilution for 10 concentration
100	10.0 ml	none
500	50.0 ml	none
1,000	50.0 ml	10.0ml vaccine plus 10.0ml diluent
2,000	50.0 ml	3.4 ml vaccine plus 10.0ml diluent
2,500	50.0 ml	2.5 ml vaccine plus 10.0ml diluent
5,000	100.0 ml	2.5 ml vaccine plus 10.0ml diluent
10,000	100.0 ml	1.2 ml vaccine plus 10.0ml diluent
25,000	100.0 ml	2.1 ml vaccine plus 50.0ml diluent

**4. Performance of the test**

**4.1 Egg inoculation**

Prior to reconstituting the vaccine, prepare and label the appropriate number of eggs for allantoic cavity inoculation.

**4.1.1 NDV inoculation**

Inoculate the  $10^{-4}$  through  $10^{-8}$  virus dilutions of the ND or ND-IB vaccine. Use 5, 9- to 11-day-old embryonated chick eggs per dilution (25 eggs total), and inoculate 0.1 ml of the appropriate dilution per embryo.

**4.1.2 IBV inoculation**

Inoculate the  $10^{-2}$  through  $10^{-6}$  dilutions of the IB or the NDV-neutralized ND-IB vaccines. Use 6, 9- to 11-day-old embryonated chick eggs (30 eggs total). Inoculate 0.1 ml of the appropriate dilution of the IB vaccine or 0.2 ml of the appropriate dilution of the NDV-neutralized ND-IB vaccine per embryo.

**4.2 Incubation**

Incubate the eggs for 7 days candling daily. Deaths occurring the first 24 hr shall be considered due to trauma and not used in calculations. At least 4 embryos per dilution must be viable at 24 hr postinoculation (PI) for a valid test.

**5. Interpretation of the test results**

**Supplemental Assay Method for the Titration of Newcastle Disease Vaccine,  
Infectious Bronchitis Vaccine, and Combination Newcastle Disease-Infectious  
Bronchitis Vaccine in Chicken Embryos**

**5.1 Controls**

Titrate a known positive reference virus with each group of titrations. The titer of the positive reference virus must be within the established range for the test results to be valid.

**5.2 NDV**

Record all deaths occurring after 24 hr as positive. On the seventh day PI, open all remaining eggs and examine the embryos. All obviously stunted embryos are considered positive.

**5.3 IBV**

Record all deaths occurring after 24 hr as positive. On the seventh day PI, open all remaining eggs and examine the embryos for IBV lesions. An embryo exhibiting 1 or more lesions is considered positive.

**1. Massachusetts type**

Check for stunting, curling, and clubbed down.

**2. Other types**

Check for stunting, curling, and clubbed down. Open the embryos and check for bile stasis (dark green liver) and kidney urates.

**5.4 Calculations**

Determine the log<sub>10</sub> EID<sub>50</sub> titer using the method of Reed and Muench. This dilution and inoculation procedure allow for a direct readout on a per dose basis. Round to 1 decimal.



Supplemental Assay Method for the Titration of Newcastle Disease Vaccine,  
Infectious Bronchitis Vaccine, and Combination Newcastle Disease-Infectious  
Bronchitis Vaccine in Chicken Embryos

## 5.5 Retests

Conduct retests as required by the Code of Federal Regulations, Title 9, Part 113.8 (b) and requirements of minimum release in firm's current Outline of Production, Part V.

## 5.6 Evaluation of test results

5.5.1 The 9CFR 113.8 (b) defines the criteria for a satisfactory/unsatisfactory serial.

5.5.2 The firm's requirements of minimum release/stability titers for each vaccine are listed in the current Outline of Production, Part V, for the specific product code.

## 6. Report of test results

Titers are reported out as EID<sub>50</sub> per bird dose.

## 7. References

1. Reed, L. J., and H. Muench. 1938. A simple method of estimating fifty percent endpoints. Am. J. Hyg. 27:493-497.
2. This document was rewritten to meet the current CVB-L QA SAM format. No significant changes were made from the previous protocol. This document supersedes the September 1, 1985, version.